



# Stereological formulas

AREA FRACTION FRACTIONATOR .....	3
CAVALIERI ESTIMATOR.....	4
COMBINED POINT INTERCEPT.....	6
CONNECTIVITY ASSAY .....	7
CYCLOIDS FOR LV.....	8
CYCLOIDS FOR SV.....	10
DISCRETE VERTICAL ROTATOR .....	12
FRACTIONATOR.....	13
ISOTROPIC FAKIR.....	17
ISOTROPIC VIRTUAL PLANES .....	18
IUR PLANES OPTICAL FRACTIONATOR .....	21
L-CYCLOID OPTICAL FRACTIONATOR .....	22
MERZ .....	23



# Stereological formulas

NUCLEATOR.....	24
OPTICAL FRACTIONATOR .....	25
OPTICAL ROTATOR.....	29
PETRIMETRICS.....	31
PHYSICAL FRACTIONATOR.....	32
PLANAR ROTATOR.....	33
POINT SAMPLED INTERCEPT .....	35
SIZE DISTRIBUTION.....	36
SPACEBALLS .....	38
SURFACE-WEIGHTED STAR VOLUME.....	40
SURFACTOR .....	41
SV-CYCLOID FRACTIONATOR .....	42
VERTICAL SPATIAL GRID .....	43
WEIBEL.....	44



## Stereological formulas

### AREA FRACTION FRACTIONATOR

<p>Estimated volume fraction (<math>\hat{V}_v</math>)</p>	$\hat{V}_v(Y, ref) = \frac{\sum_{i=1}^m P(Y)_i}{\sum_{i=1}^m P(ref)_i}$	<p><math>P(ref)</math> Points hitting reference volume  <math>Y</math> Sub-region  <math>P(Y)</math> Points hitting sub-region</p>
<p>Estimated area (<math>\hat{A}</math>)</p>	$\hat{A} = \frac{1}{asf} \cdot a(p) \cdot P(Y_i)$	<p><math>asf</math> Area sampling fraction  <math>a(p)</math> Area associated with a point</p>

### References

Howard, C. V., & Reed, M. G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (pp. 170–172). Milton Park, England: BIOS Scientific Publishers.



## Stereological formulas

### CAVALIERI ESTIMATOR

<b>Area associated with a point (<math>A_p</math>)</b>	$A_p = g^2$	$g^2$ Grid area
<b>Volume associated with a point (<math>V_p</math>)</b>	$V_p = g^2 m \bar{t}$	$m$ Section evaluation interval $\bar{t}$ Mean section cut thickness
<b>Estimated volume (<math>\hat{V}</math>)</b>	$\hat{V} = A_p m' \bar{t} \left( \sum_{i=1}^n P_i \right)$	$A_p$ Area associated with a point $m'$ Section evaluation interval $\bar{t}$ Mean section cut thickness $P_i$ Points counted on grid
<b>Estimated volume corrected for over-projection (<math>[v]</math>)</b>	$[v] = t \cdot \left( k \cdot \sum_{j=1}^g a'_j - \max(a') \right)$	$t$ Section cut thickness $k$ Correction factor $g$ Grid size $a'$ Projected area
<b>Coefficient of error (CE)</b>	$CE = \frac{\sqrt{TotalVar}}{\sum_{i=1}^n P_i}$	$TotalVar$ Total variance of the estimated volume $n$ Number of sections $P_i$ Points counted on grid  $TotalVar = s^2 + VAR_{SRS}$



# Stereological formulas

## Cavalieri Estimator (2)

<p><b>Variance of systematic random sampling</b> (<math>VAR_{SRS}</math>)</p>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	<p><math>m</math> Smoothness class of sampled function  <math>s^2</math> Variance due to noise  <math>A = \sum_{i=1}^n P_i^2, B = \sum_{i=1}^{n-1} P_i P_{i+1}, C = \sum_{i=1}^{n-2} P_i P_{i+2}</math></p> <p>With:</p> <p><math>n</math> : number of sections</p> $s^2 = 0.0724 \left(\frac{b}{\sqrt{a}}\right) \sqrt{n \sum_{i=1}^n P_i}$ <p><math>\frac{b}{\sqrt{a}}</math> Shape factor</p>
---	--	--

## References

García-Fiñana, M., Cruz-Orive, L.M., Mackay, C.E., Pakkenberg, B. & Roberts, N. (2003). [Comparison of MR imaging against physical sectioning to estimate the volume of human cerebral compartments.](#) *Neuroimage*, 18 (2), 505–516.

Gundersen, H. J. G., & Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction.](#) *Journal of Microscopy*, 147 (3), 229–263.

Howard, C. V., & Reed, M.G. (2005). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (Chapter 3). New York: Garland Science/BIOS Scientific Publishers.



## Stereological formulas

### COMBINED POINT INTERCEPT

<b>Profile area (<math>a</math>)</b>	$a = a(p) \cdot \sum P$	$a(p)$ Area associated with a point $\sum P$ Number of points
<b>Profile boundary (<math>b</math>)</b>	$b = \frac{\pi}{2} d \cdot \sum I$	$d$ Distance between points $\sum I$ Number of intersections

This method is based on the principles described in the following:

Howard, C.V., Reed, M.G. (2010). *Unbiased Stereology* (Second Edition). QTP Publications: Coleraine, UK. See equations 2.5 and 3.2

Miles, R.E., Davy, P. (1976). Precise and general conditions for the validity of a comprehensive set of stereological fundamental formulae. *Journal of Microscopy*, 107 (3), 211–226.



## Stereological formulas

### CONNECTIVITY ASSAY

<b>Euler number (<math>X_3</math>)</b>	$X_3 = I + H - B$	$I$ Total island markers $H$ Total hole markers $B$ Total bridge markers
<b>Number of alveoli (<math>N_{alv}</math>)</b>	$N_{alv} = -X_3$	$X_3$ Euler number
<b>Sum counting frame volumes (<math>V</math>)</b>	$V = h \cdot n \cdot a$	$h$ Disector height $n$ Number of disectors $a$ Area counting frame
<b>Numerical density of alveoli (<math>N_v</math>)</b>	$N_v = \frac{N_{alv}}{V}$	$N_{alv}$ Number of alveoli $V$ Sum counting frame volumes

### References

Ochs, M., Nyengaard, J.R., Jung, A., Knudsen, L., Voigt, M., Wahlers, T., Richter, J., & Gundersen, H.J.G. (2004). [The number of alveoli in the human lung](#). *American journal of respiratory and critical care medicine*, 169 (1), 120–124.



## Stereological formulas

### CYCLOIDS FOR LV

<b>Area associated with a point (<math>A_p</math>)</b>	$A_p = g^2$	$g^2$ Grid area
<b>Volume associated with a point (<math>V_p</math>)</b>	$V_p = g^2 m \bar{t}$	$g^2$ Grid area $m$ Section evaluation interval $\bar{t}$ Mean section cut thickness
<b>Length per unit volume (<math>L_V</math>)</b>	$L_V = 2 \frac{[\bar{I}_L^C]_{prj}}{\Delta}$ $L_V = \frac{2}{\Delta} \cdot \frac{(\bar{I}_C^{cyc})_{prj}}{\bar{p} \cdot \left(\frac{l}{p}\right)} = \frac{2}{\Delta} \left(\frac{p}{l}\right) \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	$[\bar{I}_L^C]_{prj}$ Number of counting frames $\Delta$ Section cut thickness $I_i$ Intercepts $P_i$ Test points $[\bar{I}_C^{cyc}]_{prj}$ Average number of intersections of projected images $\frac{p}{l}$ Test points per unit length of cycloid
<b>Estimated volume (<math>\hat{V}</math>)</b>	$\hat{V} = m \Delta \left(\frac{a}{p}\right) \sum_{i=1}^n P_i$	$m$ Sampling fractions $\Delta$ Section cut thickness $a$ Area $p$ Number of test points $P_i$ Test points
<b>Estimated length (<math>\hat{L}</math>)</b>	$\hat{L} = 2 \left(\frac{a}{l}\right) m \sum_{i=1}^n I_i$	$a$ Area $l$ Line length $m$ Sampling fractions $I_i$ Intercepts





## Stereological formulas

### Cycloids for $L_V$ (2)

<b>Coefficient of error for line length</b>	$CE(\hat{L} L) = \frac{\sqrt{VAR_{SRS}}}{\sum_{i=1}^n I_i}$	$VAR_{SRS}$ Variance of systematic random sampling $\hat{L} L$ Estimated length per length $I_i$ Intercepts
<b>Variance of systematic random sampling (<math>VAR_{SRS}</math>)</b>	$VAR_{SRS} = \frac{3g_0 - 4g_1 + g_2}{12}$ $g_k = \sum_{i=1}^{n-k} L_i L_{i+k}$	$g$ Grid size $L_i$ Line length at section $i$
<b>Coefficient of error for length density</b>	$CE(L_V) = \sqrt{\frac{n}{n-1} \left( \frac{\sum_{i=1}^n I_i^2}{\sum_{i=1}^n I_i \sum_{i=1}^n I_i} + \frac{\sum_{i=1}^n P_i^2}{\sum_{i=1}^n P_i \sum_{i=1}^n P_i} - 2 \frac{\sum_{i=1}^n I_i P_i}{\sum_{i=1}^n I_i \sum_{i=1}^n P_i} \right)}$	$I_i$ Intercepts $P_i$ Test points $n$ Number of probes

## References

Artacho-Pérula, E., Roldán-Villalobos, R. (1995). Estimation of capillary length density in skeletal muscle by unbiased stereological methods: I. Use of vertical slices of known thickness *The Anatomical Record*, 241 (3), 337-344.

Gokhale, A. M. (1990). Unbiased estimation of curve length in 3-D using vertical slices. *Journal of Microscopy*, 159 (2), 133–141.

Howard, C. V., Reed, M.G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (pp. 170–172). BIOS Scientific Publishers.



## Stereological formulas

### CYCLOIDS FOR SV

<b>Area associated with a point (<math>A_p</math>)</b>	$A_p = g^2$	$g^2$ Grid area
<b>Volume associated with a point (<math>V_p</math>)</b>	$V_p = g^2 m \bar{t}$	$g^2$ Grid area $m$ Evaluation interval $\bar{t}$ Section cut thickness
<b>Estimated surface area per unit volume (<math>est S_v</math>)</b>	$est S_v = 2 \left( \frac{2p}{l} \right) \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	$p/l$ Points per unit length of cycloid $I_i$ Intercepts with cycloids $P_i$ Point counts
<b>Estimated volume (<math>\hat{V}</math>)</b>	$\hat{V} = m \bar{t} \left( \frac{a}{p} \right) \sum_{i=1}^m P_i$	$m$ Evaluation interval $\bar{t}$ Section cut thickness $a/p$ Area associated with each point $P_i$ Point counts
<b>Estimated surface area (<math>\hat{S}</math>)</b>	$\hat{S} = 2 \left( \frac{a}{l} \right) m \bar{t} \sum_{i=1}^m I_i$	$m$ Evaluation interval $\bar{t}$ Section cut thickness $a/l$ Area per unit length $I_i$ Intercepts with cycloids



## Stereological formulas

### Cycloids for Sv (2)

<p><b>Coefficient of error for estimated surface (CE)</b></p>	$CE(\hat{S}) = \frac{\sqrt{VAR_{SRS}}}{\sum_{i=1}^n I_i}$	<p>Var<sub>SRS</sub> Variance due to systematic random sampling</p> $Var_{SRS} = \frac{3g_0 - 4g_1 + g_2}{12}$
<p><b>Coefficient of error for surface density (CE (S<sub>v</sub>))</b></p>	$CE(S_v) = \sqrt{\frac{n}{n-1} \left( \frac{\sum_{i=1}^n I_i^2}{\sum_{i=1}^n I_i \sum_{i=1}^n I_i} + \frac{\sum_{i=1}^n P_i^2}{\sum_{i=1}^n P_i \sum_{i=1}^n P_i} - 2 \frac{\sum_{i=1}^n I_i P_i}{\sum_{i=1}^n I_i \sum_{i=1}^n P_i} \right)}$	<p><i>n</i> Number of measurements  <i>I<sub>i</sub></i> Intercepts with cycloids  <i>P<sub>i</sub></i> Point counts</p>

### References

Baddeley, A. J., Gundersen, H.J.G., & Cruz-Orive, L.M. (1998) Estimation of surface area from vertical sections. *Journal of Microscopy*, 142 (3), 259–276.

Howard, C. V., Reed, M.G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy*(pp.170–172). BIOS Scientific Publishers.



## Stereological formulas

### DISCRETE VERTICAL ROTATOR

<b>Estimated volume</b> ( <i>Est v</i> )	$est\ v = \frac{\pi}{n} \cdot a_p \cdot \sum_{i=1}^n P_i D_i$	<i>n</i> Number of centriolar sections <i>a<sub>p</sub></i> Area associated with each point <i>P<sub>i</sub></i> Number of points in each class <i>D<sub>i</sub></i> Distance of class from central axis
---	---	---

### References

Mironov, A. A. (1998). Estimation of subcellular organelle volume from ultrathin sections through centrioles with a discretized version of the vertical rotator. *Journal of microscopy*, 192(1), 29-36.



## Stereological formulas

### FRACTIONATOR

<b>Estimate of total number of particles (<math>N</math>)</b>	$N = \sum Q^- \cdot \frac{1}{asf} \cdot \frac{1}{ssf}$	$Q^-$ Particles counted $asf$ Area sampling fraction $ssf$ Section sampling fraction
<b>Variance due to systematic random sampling – Gundersen (<math>VAR_{SRS}</math>)</b>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-$ $s^2$ Variance due to noise
<b>Variance due to noise – Gundersen (<math>s^2</math>)</b>	$s^2 = \sum_{i=1}^n Q^-$	$Q^-$ Particles counted $n$ Number of sections used
<b>Total variance – Gundersen (<math>TotalVar</math>)</b>	$TotalVar = s^2 + VAR_{SRS}$	$VAR_{SRS}$ Variance due to SRS $s^2$ Variance due to noise
<b>Coefficient of error – Gundersen (<math>CE</math>)</b>	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$ Total variance $s^2$ Variance due to noise
<b>Number-weighted mean section cut thickness (<math>\bar{t}_{Q^-}</math>)</b>	$\bar{t}_{Q^-} = \frac{\sum_{i=1}^m t_i Q_i^-}{\sum_{i=1}^m Q_i^-}$	$m$ Number of sections $t_i$ Section thickness at site $i$ $Q_i$ Particles counted



## Stereological formulas

### Fractionator (2)

<b>Coefficient of error – Scheaffer (CE)</b>	$CE = \frac{\sqrt{s^2 \left( \frac{1}{f} - \frac{1}{F} \right)}}{\bar{Q}}$	<p><math>f</math> Number of counting frames  <math>F</math> Total possible sampling sites  <math>s^2</math> Estimated variance  <math>\bar{Q}</math> Average particles counted</p>
<b>Average number of particles – Scheaffer (<math>\bar{Q}</math>)</b>	$\bar{Q} = \frac{\sum_{i=1}^f Q_i}{f}$	<p><math>Q_i</math> Particles counted  <math>f</math> Number of counting frames</p>
<b>Estimated variance - Scheaffer (<math>s^2</math>)</b>	$s^2 = \frac{\sum_{i=1}^f (Q_i - \bar{Q})^2}{f - 1}$	<p><math>f</math> Number of counting frames  <math>Q_i</math> Particles counted  <math>\bar{Q}</math> Average particles counted</p>
<b>Estimated variance of estimated cell population - Scheaffer</b>	$\frac{C_{fp} F^2 s^2}{f}$	<p><math>C_{fp}</math> Finite population correction  <math>s^2</math> Estimated variance  <math>f</math> Number of counting frames  <math>F</math> Total possible sampling sites</p>
<b>Estimated variance of mean cell count - Scheaffer</b>	$\frac{C_{fp} s^2}{f}$	<p><math>C_{fp}</math> Finite population correction  <math>s^2</math> Estimated variance  <math>f</math> Number of counting frames</p>



## Stereological formulas

### Fractionator (3)

<p><b>Estimated mean coefficient of error – Cruz-Orive (<i>est Mean CE</i>)</b></p>	$est\ Mean\ CE\ (est\ N) = \left[ \frac{1}{3n} \cdot \sum_{i=1}^n \left( \frac{Q_{1i} - Q_{2i}}{Q_{1i} + Q_{2i}} \right)^2 \right]^{1/2}$	<p><math>Q_{1i}</math> Counts in sub-sample 1  <math>Q_{2i}</math> Counts in sub-sample 2  <math>n</math> Size of sub-sample</p>
<p><b>Predicted coefficient of error for estimated population – Schmitz-Hof (<math>CE_{pred}</math>)</b></p>	$CE_{pred}(n_F) = \sqrt{\frac{Var(Q_r^-)}{R \cdot (Q_r^-)^2}}$ $CE_{pred}(n_F) = \frac{1}{\sqrt{\sum_{r=1}^R Q_r^-}} = \frac{1}{\sqrt{\sum_{s=1}^S Q_s^-}}$	<p><math>R</math> Number of counting spaces  <math>S</math> Number of sections  <math>Q_r^-</math> Counts in the "r"-th counting space  <math>Q_s^-</math> Counts in the "s"-th section</p>

### References

- Geiser, M., Cruz-Orive, L.M., Hof, V.I., & Gehr, P. (1990) Assessment of particle retention and clearance in the intrapulmonary conducting airways of hamster lungs with the fractionator. *Journal of Microscopy*, 160 (1), 75–88.
- Glaser, E. M., Wilson, P.D. (1998). The coefficient of error of optical fractionator population size estimates: a computer simulation comparing three estimators. *Journal of Microscopy*, 192 (2), 163–171.
- Gundersen, H.J.G., Vedel Jensen, E.B., Kieu, K., & Nielsen, J. (1999). The efficiency of systematic sampling in stereology—reconsidered. *Journal of Microscopy*, 193 (3), 199–211.
- Gundersen, H. J. G., Jensen, E.B. (1987). The efficiency of systematic sampling in stereology and its prediction. *Journal of Microscopy*, 147 (3), 229–263.
- Howard, V., Reed, M. (2005). *Unbiased stereology: three-dimensional measurement in microscopy* (vol. 4, chapter 12). Garland Science/Bios Scientific Publishers.



## Stereological formulas

### Fractionator (4)

Scheaffer, R.L., Ott, L., & Mendenhall, W. (1996). *Elementary survey sampling* (chapter 7). Boston: PWS-Kent.

Schmitz, C., Hof, P.R. (2000). Recommendations for straightforward and rigorous methods of counting neurons based on a computer simulation approach. *Journal of Chemical Neuroanatomy*, 20 (1), 93–114.

West, M. J., Slomianka, L., & Gundersen, H.J.G. (1991). Unbiased stereological estimation of the total number of neurons in the subdivisions of the rat hippocampus using the optical fractionator. *The Anatomical Record*, 231 (4), 482–497.





## Stereological formulas

### ISOTROPIC FAKIR

<b>Estimated total surface area</b>	$estS = 2 \frac{1}{n} \cdot \sum_{i=1}^n \frac{v}{l_i} \cdot I_i$	$n$ Number of line sets (always set to 3) $\frac{v}{l_i}$ Inverse of the probe per unit volume $I_i$ Intercepts with test lines
-------------------------------------	---	---

### References

Kubínová, L., Janacek, J. (1998). Estimating surface area by the isotropic fakir method from thick slices cut in an arbitrary direction. *Journal of Microscopy*, 191 (2), 201–211.



## Stereological formulas

### ISOTROPIC VIRTUAL PLANES

<b>Length per unit volume</b>	$L_v = \frac{2p(box)}{a(plane)} \cdot \frac{\sum Q}{\sum p(ref)}$	<p><math>p(box)</math> Number of corners considered</p> <p><math>a(plane)</math> Exact sampling area</p> <p><math>p(ref)</math> Number of corners in region</p> <p><math>\sum Q</math> Total number of transects</p>
<b>Estimated total length</b>	$L = \frac{1}{ssf} \cdot \frac{1}{asf} \cdot \frac{1}{hsf} \cdot \frac{1}{psd} \cdot 2 \sum Q$ <p>Or <math>L = \frac{1}{ssf} \cdot \frac{dx \cdot dy}{a(box)} \cdot \frac{\bar{t}}{h(box)} \cdot d \cdot 2 \sum Q</math></p> $asf = \frac{a(box)}{dx \cdot dy}$ $hsf = \frac{h(box)}{\bar{t}}$ $psd = \frac{E[a(plane)]}{v(box)} = \frac{1}{d}$	<p><math>ssf</math> Section sampling fraction</p> <p><math>asf</math> Area sampling fraction</p> <p><math>hsf</math> Height sampling fraction</p> <p><math>psd</math> Probe sampling density</p> <p><math>\sum Q</math> Total number of transects</p> <p><math>a(box)</math> Area of sampling box</p> <p><math>h(box)</math> Depth of sampling box</p> <p><math>d</math> Sampling plane separation</p> <p><math>dx, dy</math> Distances in XY</p> <p><math>\bar{t}</math> Average section thickness</p> <p><math>a(plane)</math> Sampling plane area</p> <p><math>E</math> Expected value</p> <p><math>v(box)</math> Volume of sampling box</p>
<b>Total plane area</b>	$A = \sum_{i=1}^l \sum_{j=1}^s A_{i,j}$	<p><math>l</math> Number of layouts</p> <p><math>s</math> Number of sampling sites</p> <p><math>A_{ij}</math> Plane area inside of each sampling box for each layout</p>



# Stereological formulas

## Isotropic Virtual Planes (2)

<p><b>Coefficient of error (Gundersen)</b></p>	$CE = \frac{\sqrt{\frac{3(A - \text{Poisson noise})}{12} + \text{Poisson noise}}}{\sum Q}, m = 0$ $CE = \frac{\sqrt{\frac{3(A - \text{Poisson noise}) - 4B + C}{240} + \text{Poisson noise}}}{\sum Q}, m = 1$	<p>Poisson noise is <math>\sum Q_i</math></p> $A = \sum (Q_i \cdot Q_i)$ $B = \sum (Q_i \cdot Q_{i+1})$ $C = \sum (Q_i \cdot Q_{i+2})$ <p><math>A, B, C</math> Covariogram values  <math>Q_i</math> Particles counted</p>
<p><b>Plane areas</b></p>	$\vec{V} = (A, B, C)$ <p>planes: <math>Ax + By + Cz + D_i = 0, \quad D_i = D + d \cdot i</math></p> $\text{box: } \left\{ \left\{ (x, y, z) \left  \begin{array}{l} x_0 \leq x \leq x_0 + b_x \\ y_0 \leq y \leq y_0 + b_y \\ z_0 \leq z \leq z_0 + b_z \end{array} \right. \right\} \right\}$ <p>area: <math>(\text{box} \cap \text{planes})</math></p>	<p><math>A, B, C, D</math> Given constants  <math>d</math> Distance between planes  <math>i</math> Integer  <math>x_0, y_0, z_0</math> Vertex of a sampling box  <math>b_x, b_y, b_z</math> Dimensions of a sampling box</p>



## Stereological formulas

### Isotropic Virtual Planes (3)

<b>Average number of counts</b>	$\bar{Q} = \frac{\sum_{i=1}^p \sum_{j=1}^{l_j} Q_{ij}}{\sum_{j=1}^p l_j}$	$p$ Number of probes $l_j$ Number of layouts in each probe $Q_{ij}$ Number of counts in each probe and layout
<b>Total corners of sampling boxes inside the region of interest</b>	$C = \sum_{i=1}^p C_i$	$p$ Number of probes $C_i$ Number of sampling boxes inside region of interest

### References

Larsen, J. O., Gundersen, H.J.G., & Nielsen, J. (1998). Global spatial sampling with isotropic virtual planes: estimators of length density and total length in thick, arbitrarily orientated sections. *Journal of Microscopy*, 191, 238–248.



## Stereological formulas

### IUR PLANES OPTICAL FRACTIONATOR

<p><b>Estimated length</b></p>	$est L = 2 \cdot \frac{a}{l} \sum I \cdot \frac{1}{ssf} \cdot \frac{1}{asf} \cdot \frac{t}{h}$	<p><math>a/l</math> Area per unit length of test line  <math>\sum I</math> Number of intersections  <math>ssf</math> Section sampling fraction  <math>asf</math> Area sampling fraction  <math>t</math> Section cut thickness  <math>h</math> Height of counting frame</p>
<p><b>Area sampling fraction</b></p>	$asf = \frac{area(Frame)}{area(x, y step)} = \frac{x_{CF} \cdot y_{CF}}{x_{step} \cdot y_{step}}$	<p><math>x_{CF}, y_{CF}</math> XY dimensions of counting frame  <math>x_{step}, y_{step}</math> Dimensions of grid  <math>area(Frame)</math> Area of counting frame  <math>area(x, y step)</math> Area of grid</p>



## Stereological formulas

### L-CYCLOID OPTICAL FRACTIONATOR

<p><b>Estimated length of lineal structure</b></p>	$est L = 2 \cdot \frac{a}{l} \sum I \cdot \frac{1}{ssf} \cdot \frac{1}{asf} \cdot \frac{t}{h}$	<p><math>a/l</math> Area per unit cycloid length  <math>\sum I</math> Number of intercepts  <math>ssf</math> Section sampling fraction  <math>asf</math> Area sampling fraction  <math>t</math> Section cut thickness  <math>h</math> Height of counting frame</p>
<p><b>Area sampling fraction</b></p>	$asf = \frac{area(Frame)}{area(x, y step)} = \frac{x_{CF} \cdot y_{CF}}{x_{step} \cdot y_{step}}$	<p><math>x_{CF}, y_{CF}</math> XY dimensions of counting frame  <math>x_{step}, y_{step}</math> Dimensions of grid  <math>area(Frame)</math> Area of counting frame  <math>area(x, y step)</math> Area of grid</p>

### References

Stocks, E. A., McArthur, J.C., Griffen, J.W., & Mouton, P.R. (1996). An unbiased method for estimation of total epidermal nerve fiber length. *Journal of Neurocytology*, 25 (1), 637–644.



## Stereological formulas

### MERZ

<b>Length of semi-circle (<math>L</math>)</b>	$L = \frac{1}{2}\pi d$	$d$ Circle diameter
<b>Surface area per unit volume (<math>S_v</math>)</b>	$S_v = \frac{2\sum I}{\frac{l}{p}\sum P}$	$I$ Number of intercepts $l/p$ Length of semi-circle per point $P$ Number of points

### References

Howard, C. V., Reed, M. G. (2010). *Unbiased stereology*. Liverpool, UK: QTP Publications. {See equation 6.4}

Weibel, E.R. (1979). *Stereological Methods. Vol. 1: Practical methods for biological morphometry*. London, UK: Academic Press.



## Stereological formulas

### NUCLEATOR

<b>Area estimate</b>	$a = \pi \bar{l}^2$	$l$ Length of rays
<b>Volume estimate</b>	$\bar{v}_N = \frac{4\pi}{3} \bar{l}_n^3$	$l$ Length of rays
<b>Estimated coefficient of error</b>	$est\ CV(R) = \frac{\sqrt{\frac{1}{n-1} \sum_{i=1}^n (R_i - \bar{R})^2}}{\bar{R}}$	$n$ Number of nucleator estimates $R_i$ Area/volume estimate for each sampling site
<b>Average area/volume estimate</b>	$\bar{R} = \frac{1}{n} \sum_{i=1}^n R_i$	$n$ Number of nucleator estimates $R_i$ Area/volume estimate for each sampling site
<b>Relative efficiency</b>	$CE_n(R) = \frac{CV(R)}{\sqrt{n}}$	$n$ Number of nucleator estimates $CV(R)$ Estimated coefficient of variation
<b>Geometric mean of area/volume estimate</b>	$e^{\left(\frac{1}{n} \sum_{i=1}^n \ln R_i\right)}$	$n$ Number of nucleator estimates $R_i$ Area/volume estimate for each sampling site

### References

Gundersen, H.J.G. (1988). The nucleator. *Journal of Microscopy*, 151 (1), 3–21.





## Stereological formulas

### OPTICAL FRACTIONATOR

<b>Estimate of total number of particles (<math>N</math>)</b>	$N = \sum Q^- \cdot \frac{t}{h} \cdot \frac{1}{asf} \cdot \frac{1}{ssf}$	$Q^-$ Particles counted $t$ Section mounted thickness $h$ Counting frame height $asf$ Area sampling fraction $ssf$ Section sampling fraction
<b>Variance due to systematic random sampling – Gundersen (<math>VAR_{SRS}</math>)</b>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-$ $s^2$ Variance due to noise
<b>Variance due to noise - Gundersen (<math>s^2</math>)</b>	$s^2 = \sum_{i=1}^n Q^-$	$Q^-$ Particles counted $n$ Number of sections used
<b>Total variance – Gundersen (<math>TotalVar</math>)</b>	$TotalVar = s^2 + VAR_{SRS}$	$VAR_{SRS}$ Variance due to SRS $s^2$ Variance due to noise
<b>Coefficient of error – Gundersen (<math>CE</math>)</b>	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$ Total variance $s^2$ Variance due to noise
<b>Number-weighted mean section cut thickness (<math>t_{Q^-}</math>)</b>	$\overline{t_{Q^-}} = \frac{\sum_{i=1}^m t_i Q_i^-}{\sum_{i=1}^m Q_i^-}$	$m$ Number of sections $t_i$ Section thickness at site $i$ $Q_i$ Particles counted



## Stereological formulas

### Optical Fractionator (2)

<b>Coefficient of error – Scheaffer (CE)</b>	$CE = \frac{\sqrt{s^2 \left( \frac{1}{f} - \frac{1}{F} \right)}}{\bar{Q}}$	<i>f</i> Number of counting frames <i>F</i> Total possible sampling sites <i>s</i> <sup>2</sup> Estimated variance $\bar{Q}$ Average particles counted
<b>Average number of particles – Scheaffer (<math>\bar{Q}</math>)</b>	$\bar{Q} = \frac{\sum_{i=1}^f Q_i}{f}$	$Q_i$ Particles counted <i>f</i> Number of counting frames
<b>Estimated variance - Scheaffer (<i>s</i><sup>2</sup>)</b>	$s^2 = \frac{\sum_{i=1}^f (Q_i - \bar{Q})^2}{f - 1}$	<i>f</i> Number of counting frames $Q_i$ Particles counted $\bar{Q}$ Average particles counted
<b>Estimated variance of estimated cell population - Scheaffer</b>	$\frac{C_{fp} F^2 s^2}{f}$	$C_{fp}$ Finite population correction $s^2$ Estimated variance <i>f</i> Number of counting frames <i>F</i> Total possible sampling sites
<b>Estimated variance of mean cell count - Scheaffer</b>	$\frac{C_{fp} s^2}{f}$	$C_{fp}$ Finite population correction $s^2$ Estimated variance <i>f</i> Number of counting frames



## Stereological formulas

### Optical Fractionator (3)

<p><b>Estimated mean coefficient of error – Cruz-Orive (<i>est Mean CE</i>)</b></p>	$est\ Mean\ CE\ (est\ N) = \left[ \frac{1}{3n} \cdot \sum_{i=1}^n \left( \frac{Q_{1i} - Q_{2i}}{Q_{1i} + Q_{2i}} \right)^2 \right]^{1/2}$	<p><math>Q_{1i}</math> Counts in sub-sample 1  <math>Q_{2i}</math> Counts in sub-sample 2  <math>n</math> Size of sub-sample</p>
<p><b>Predicted coefficient of error for estimated population – Schmitz-Hof (<math>CE_{pred}</math>)</b></p>	$CE_{pred}(n_F) = \sqrt{\frac{Var(Q_r^-)}{R \cdot (Q_r^-)^2}}$ $CE_{pred}(n_F) = \frac{1}{\sqrt{\sum_{r=1}^R Q_r^-}} = \frac{1}{\sqrt{\sum_{s=1}^S Q_s^-}}$	<p><math>R</math> Number of counting spaces  <math>S</math> Number of sections  <math>Q_r^-</math> Counts in the "r"-th counting space  <math>Q_s^-</math> Counts in the "s"-th section</p>

### References

Geiser, M., Cruz-Orive, L.M., Hof, V.I., & Gehr, P. (1990) Assessment of particle retention and clearance in the intrapulmonary conducting airways of hamster lungs with the fractionator. *Journal of Microscopy*, 160 (1), 75–88.

Glaser, E. M., Wilson, P.D. (1998). The coefficient of error of optical fractionator population size estimates: a computer simulation comparing three estimators. *Journal of Microscopy*, 192 (2), 163–171.

Gundersen, H.J.G., Vedel Jensen, E.B., Kieu, K., & Nielsen, J. (1999). The efficiency of systematic sampling in stereology—reconsidered. *Journal of Microscopy*, 193 (3), 199–211.

Gundersen, H. J. G., Jensen, E.B. (1987). The efficiency of systematic sampling in stereology and its prediction. *Journal of Microscopy*, 147 (3), 229–263.

Howard, V., Reed, M. (2005). *Unbiased stereology: three-dimensional measurement in microscopy* (vol. 4, chapter 12). Garland Science/Bios Scientific Publishers.



## Stereological formulas

### Optical Fractionator (4)

Scheaffer, R.L., Ott, L., & Mendenhall, W. (1996). *Elementary survey sampling* (chapter 7). Boston: PWS-Kent.

Schmitz, C., Hof, P.R. (2000). Recommendations for straightforward and rigorous methods of counting neurons based on a computer simulation approach. *Journal of Chemical Neuroanatomy*, 20 (1), 93–114.

West, M. J., Slomianka, L., & Gundersen, H.J.G. (1991). Unbiased stereological estimation of the total number of neurons in the subdivisions of the rat hippocampus using the optical fractionator. *The Anatomical Record*, 231 (4), 482–497.



# Stereological formulas

## OPTICAL ROTATOR

<p><b>Volume of particle</b></p>	$\hat{v} = a \sum_i^{+/-} g(P_i)$	<p><math>a</math> Reciprocal line density</p> <p><math>a=k.h</math>  <math>k</math> Length of slice  <math>h</math> Systematic spacing</p>
<p><b>For vertical slabs and lines parallel to vertical axis</b></p>	$g(P) = d_1, \text{ if } d_2 < t$ $g(P) = \frac{\frac{\pi}{2} d_1}{\arcsin\left(\frac{t}{d_2}\right)}, \text{ if } t \leq d_2$	<p><math>d_1</math> Distance along test line  <math>d_2</math> Distance from origin to test line  <math>t</math> <math>\frac{1}{2}</math> thickness of optical slice</p>
<p><b>For vertical slabs and lines perpendicular to vertical axis</b></p>	$g(P) = d_1, \text{ if } \sqrt{d_1^2 + z^2} < t$ $g(P) = f\left(\sqrt{t^2 - z^2}\right), \text{ if }  z  < t \leq \sqrt{d_1^2 + z^2}$ $g(P) = f(0), \text{ if } t \leq  z $ $f(x) = x + \frac{\pi}{2} \int_x^{d_1} \frac{1}{\arcsin\left(\frac{t}{\sqrt{u^2 + z^2}}\right)} du$	<p><math>d_1</math> Distance along test line  <math>t</math> <math>\frac{1}{2}</math> thickness of optical slice  <math>z</math> Distance in <math>z</math> from intercept to origin</p>



## Stereological formulas


### Optical Rotator (2)

<p><b>For isotropic slabs</b></p>		$g(P) = d_1, \quad \text{if } d_3 < t$ $g(P) = \frac{1}{2t} [h(t, d_2) + k(t, d_1, d_2, d_3)], \text{ if } d_2 < t \leq d_3$ $h(t, d) = t^2 \sqrt{1 - \frac{d^2}{t^2}}$ $k(t, d_1, d_2, d_3) = d_1 d_3 + d_2^2 \log \left( \frac{d_1 + d_3}{t + \sqrt{t^2 - d_2^2}} \right)$	<p><math>d_1</math> Distance along test line  <math>d_2</math> Distance from origin to test line  <math>d_3</math> Distance from intercept to origin  <math>t</math> <math>\frac{1}{2}</math> thickness of optical slice</p>
<p><b>Estimated surface area</b></p>		$\hat{S} = a \sum_j l_j g(l_j)$ $g(l) = 2, \quad \text{if } d_2 < t$ $g(l) = \pi \cdot \frac{1}{\arcsin\left(\frac{t}{d_2}\right)}, \quad \text{if } t \leq d_2$	<p><math>a</math> Reciprocal line density  <math>l_j</math> Number of intersections between grid line and cell boundary  <math>d_2</math> Distance from origin to test line  <math>t</math> <math>\frac{1}{2}</math> thickness of optical slice</p>

### References

Tandrup, T., Gundersen, H.J.G., & Vedel Jensen, E.B. (1997). The optical rotator *Journal of microscopy*, 186 (2), 108–120.

## PETRIMETRICS

<p><b>Total length (<math>\hat{L}</math>)</b></p>	$\hat{L} = \frac{\pi}{2} \cdot \frac{a}{l} \cdot \frac{1}{asf} \cdot \sum I$ $\hat{L} = d \cdot \frac{1}{asf} \cdot \sum I$	<p><math>a/l = 2d/\pi</math> Grid constant (<math>2d/\pi</math> units or ratio of area to length of semi-circle probe)</p> <p><math>asf</math> Area fraction (ratio of area of counting frame to grid-step)</p> <p><math>I</math> Number of intersections counted</p> <p><math>d = 2 \cdot \text{Merz-radius}</math> where the Merz-radius refers to the radius of the semi-circle used to probe.</p> 
---	---	--

## References

Howard, C. V., & Reed, M. G. (2005). *Unbiased stereology*. New York: Garland Science (prev. BIOS Scientific Publishers).



## Stereological formulas

### PHYSICAL FRACTIONATOR

<b>Total number of particles (<math>N</math>)</b>	$N = \sum Q^- \cdot \frac{1}{asf} \cdot \frac{1}{ssf}$	$Q^-$ Particles counted $asf$ Area sampling fraction $ssf$ Section sampling fraction
<b>Variance due to noise (<math>s^2</math>)</b>	$s^2 = \sum_{i=1}^n Q^-$	$Q^-$ Particles counted $n$ Number of sections used
<b>Variance due to systematic random sampling (<math>VAR_{SRS}</math>)</b>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-$ $s^2$ Variance due to noise
<b>Total variance (<math>TotalVar</math>)</b>	$TotalVar = s^2 + VAR_{SRS}$	$VAR_{SRS}$ Variance due to SRS $s^2$ Variance due to noise
<b>Coefficient of error (<math>CE</math>)</b>	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$ Total variance $s^2$ Variance due to noise

### References

Gundersen, Hans-Jørgen G. "[Stereology of arbitrary particles\\*](#)." Journal of Microscopy 143, no. 1 (1986): 3-45.

Sterio, D. C. "[The unbiased estimation of number and sizes of arbitrary particles using the disector.](#)" Journal of Microscopy 134, no. 2 (1984): 127-136.





## Stereological formulas

### PLANAR ROTATOR

<b>Volume for isotropic planar rotator</b>	$V = 2t \sum_i g_i$	<i>t</i> Separation between test lines <i>g<sub>i</sub></i> Isotropic planar rotator function
<b>Volume for vertical planar rotator</b>	$V = \pi t \sum_i l_i^2$	<i>t</i> Separation between test lines <i>l<sub>i</sub></i> Intercept length along a test line
<b>Isotropic planar rotator function</b>	$g_i(l) = l \sqrt{l^2 + a_i^2} + a_i^2 \ln \left[ \frac{l}{a_i} + \sqrt{\left(\frac{l}{a_i}\right)^2 + 1} \right]$ $g_{i+} = \sum_{j \text{ even}} g_i(l_{i j+}) - \sum_{j \text{ odd}} g_i(l_{i j+})$ $g_{i-} = \sum_{j \text{ even}} g_i(l_{i j-}) - \sum_{j \text{ odd}} g_i(l_{i j-})$ $g_i = \frac{1}{2}(g_{i+} + g_{i-})$	<i>l</i> Intercept length along a test line <i>a<sub>i</sub></i> Distance from origin to test line <i>j</i> Number of grid lines <i>l<sub>ij</sub></i> Number of intersections between the <i>j</i> -th grid line and the cell boundary



## Stereological formulas

### Planar Rotator (2)

<b>Isotropic planar rotator function (cont'd)</b>	$l_{i+}^2 = \sum_{j \text{ even}} l_{ij+}^2 - \sum_{j \text{ odd}} l_{ij+}^2$ $l_{i-}^2 = \sum_{j \text{ even}} l_{ij-}^2 - \sum_{j \text{ odd}} l_{ij-}^2$ $l_i^2 = \frac{1}{2}(l_{i+}^2 + l_{i-}^2)$	<p><math>l</math> Intercept length along a test line <math>a_i</math> Distance from origin to test line <math>j</math> Number of grid lines <math>l_{ij}</math> Number of intersections between the <math>j</math>-th grid line and the cell boundary</p>
---	--	---

### References

Jensen Vedel, E.B., Gundersen, H.J.G. (1993). The rotator *Journal of Microscopy*, 170 (1), 35–44.



## Stereological formulas

### POINT SAMPLED INTERCEPT

<b>Volume based on intercept length (<math>\hat{V}_V</math>)</b>	$\hat{V}_V = \frac{\pi}{3} \bar{l}_0^3 = \frac{\pi}{3n} \sum_{i=1}^n l_{0,i}^3$	$n$ Number of intercepts $l$ Intercept length
<b>Volume-weighted mean volume (<math>\bar{v}_V</math>)</b>	$\bar{v}_V = \frac{\sum_{i=1}^n \bar{l}_0^3}{n} \cdot \frac{\pi}{3}$	$n$ Number of intercepts $l$ Intercept length
<b>Coefficient of error (CE)</b>	$CE(\bar{l}_0^3) = \frac{\sqrt{\frac{\sum_{i=1}^n (\bar{l}_0^3)^2}{n}}}{\left(\frac{\sum_{i=1}^n \bar{l}_0^3}{n}\right)^2} - \frac{1}{n}$	$n$ Number of intercepts $l$ Intercept length
<b>Coefficient of variance (CV)</b>	$CE(\bar{l}_0^3) = CE(\bar{v}_V) \cdot \sqrt{n}$	$n$ Number of intercepts $l$ Intercept length $\bar{v}_V$ Volume-weighted mean volume
<b>Variance (<math>Variance_V</math>)</b>	$Variance_V(v) = \left[ \frac{\pi}{3} \cdot SD(\bar{l}_0^3) \right]^2 = [CV(\bar{l}_0^3) \cdot \bar{v}_V]^2$	$L$ Intercept length $\bar{v}_V$ Volume-weighted mean volume CV Coefficient of variance

### References

Gundersen, H.J.G., Jensen. E.B. (1985). Stereological Estimation of the Volume-Weighted Mean Volume of Arbitrary Particles Observed on Random Sections. *Journal of Microscopy*, 138, 127–142.

Sørensen, F.B. (1991). Stereological estimation of the mean and variance of nuclear volume from vertical sections. *Journal of Microscopy*, 162 (2), 203–229.



## Stereological formulas

### SIZE DISTRIBUTION

<b>Volume-weighted mean particle volume</b>	$\bar{v}_V = \bar{v}_N \cdot [1 + CV_N^2(v)]$	$\bar{v}_N$ Number-weighted mean volume $CV_N(v)$ Coefficient of variation
<b>Number-weighted mean volume</b>	$\bar{v}_N = \frac{\sum S}{\sum R}$ $S = Q \cdot R$	$R$ Number of contours $Q$ Number of points per contour
<b>Variance</b>	$Var_N(v) = \frac{\left[ \sum T - \frac{(\sum S)^2}{\sum R} \right] \cdot v(p)^2}{\sum R - 1}$ $T = Q^2 \cdot R$	$R$ Number of contours $Q$ Number of points per contour $v(p)$ Volume associated with a point
<b>Standard deviation</b>	$SD_N(v) = \sqrt{Var_N(v)}$ $SD_N(v) = \sqrt{\bar{v}_N \cdot (\bar{v}_V - \bar{v}_N)}$	$\bar{v}_N$ Number-weighted mean volume $\bar{v}_V$ Volume-weighted particle volume $Var_N$ Variance
<b>Coefficient of variation</b>	$CV_N(v) = \frac{SD_N(v)}{\bar{v}_N}$ $CV_N(v) = \sqrt{\frac{\bar{v}_V - \bar{v}_N}{\bar{v}_N}}$	$\bar{v}_N$ Number-weighted mean volume $\bar{v}_V$ Volume-weighted particle volume $SD_N$ Standard deviation



## Stereological formulas

Size distribution (2)

<b>Coefficient of error</b>	$CE_N(v) = \frac{CV_N(v)}{\sqrt{R}}$	$CV_N$ Coefficient of variation $R$ Number of contours
-----------------------------	--------------------------------------	---

### References

Sørensen, F.B. (1991). Stereological estimation of the mean and variance of nuclear volume from vertical sections. *Journal of Microscopy*, 162 (2), 203–229.



# Stereological formulas

## SPACEBALLS

<p><b>Length estimate</b></p>	$L = 2 \cdot \left( \sum_{i=1}^n Q_i \right) \cdot \frac{v}{a} \cdot \frac{1}{ssf}$ <p><i>This equation does not include the terms F2 (area-fraction) and F3 (thickness-fraction) used by Mouton et al. (equation 2, 2002), but includes that information in v (volume sampled).</i></p>	<p><math>n</math> Number of sections used  <math>Q_i</math> Intersection counted  <math>v</math> Volume (grid X * grid Y * section thickness)  <math>a</math> Surface area of the sphere  <math>ssf</math> Section sampling fraction</p>
<p><b>Variance due to noise</b></p>	$s^2 = \sum_{i=1}^n Q_i$	<p><math>Q_i</math> Intersection counted</p>
<p><b>Variance due to systematic random sampling</b></p>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	<p><math>A = \sum_{i=1}^n (Q_i^-)^2</math>  <math>B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-</math>  <math>C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-</math></p> <p><math>s^2</math> Variance due to noise  <math>m</math> Smoothness class of sampled function</p>
<p><b>Total variance</b></p>	$TotalVar = s^2 + VAR_{SRS}$	<p><math>VAR_{SRS}</math> Variance due to SRS  <math>s^2</math> Variance due to noise</p>



## Stereological formulas

### Spaceballs (2)

<b>Coefficient of error</b>	$CE = \frac{\sqrt{TotalVar}}{s^2}$	<i>TotalVar</i> Total variance <i>s</i> <sup>2</sup> Variance due to noise
-----------------------------	------------------------------------	---

### References

Mouton, P. R., Gokhale, A.M., Ward, N.L., & West, M.J. (2002). Stereological length estimation using spherical probes. *Journal of Microscopy*, 206 (1), 54–64.



SURFACE-WEIGHTED STAR VOLUME

<p><b>Surface-weighted star volume (<math>\hat{v}_s^*</math>)</b></p>	$\hat{v}_s^* = \frac{2\pi}{3} \cdot \bar{l}_1^3$ $\hat{v}_s^* = \frac{2\pi}{3} \cdot \frac{\sum_{i=1}^n \sum_{j=1}^{m_i} l_{1,(i,j)}^3}{\sum_{i=1}^n m_i}$	<p><math>n</math> Number of probes  <math>l</math> Intercept length  <math>m_i</math> Number of intercepts</p>
<p><b>Sum of cubed intercepts in probe (<math>y_i</math>)</b></p>	$y_i = \sum_{j=1}^{m_i} l_{1,(i,j)}^3$	<p><math>m_i</math> Number of intercepts  <math>l</math> Intercept length</p>
<p><b>Coefficient of error (CE)</b></p>	$CE[\hat{v}_s^*] = \left[ \frac{n}{n-1} \left\{ \frac{\sum y_i^2}{\sum y_i \sum y_i} + \frac{\sum m_i^2}{\sum m_i \sum m_i} - 2 \cdot \frac{\sum m_i y_i}{\sum y_i \sum m_i} \right\} \right]^{1/2}$	<p><math>n</math> Number of probes  <math>y_i</math> Sum of cubed intercepts in probe  <math>m_i</math> Number of intercepts</p>

References

Reed, M. G., Howard, C.V. (1998). Surface-weighted star volume: concept and estimation. *Journal of Microscopy*, 190 (3), 350–356.





# Stereological formulas

## SURFACTOR

<b>Surface area for single-ray designs</b>	$\hat{S} = 4\pi l_0^2 + c(\beta)$	$l$ Length of intercept $\beta$ Angle between test line and surface $c(\beta)$ Function of the planar angle
<b>Surface area for multi-ray designs</b>	$\hat{S} = 2\pi \sum_{j=1}^{2r} l_j^2 \cdot c(\beta)$	$l$ Length of intercept $\beta$ Angle between test line and surface $c(\beta)$ Function of the planar angle $r$ Number of test lines
<b>Function of the planar angle</b>	$c(\beta) = 1 + \left[ \frac{1}{2} \cot \beta \right] \cdot \left[ \frac{\pi}{2} - \sin^{-1} \frac{1 - \cot^2 \beta}{1 + \cot^2 \beta} \right]$	$\beta$ Angle between test line and surface

### References

Jensen, E.B., Gundersen, H.J.G. (1987). Stereological estimation of surface area of arbitrary particles. *Acta Stereologica*, 6 (3).



## Stereological formulas

### SV-CYCLOID FRACTIONATOR

<b>Estimated surface area per unit volume</b>	$S_V = 2 \left( \frac{p}{l} \right) \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	$p/l$ Ratio of test points to curve length $n$ Number of micrographs $\sum I_i$ Total intercept points on curve $\sum P_i$ Total test points
---	--	---

#### References

Baddeley, A. J., Gundersen, H.J.G., & Cruz-Orive, L.M. (1986). Estimation of surface area from vertical sections. *Journal of Microscopy*, 142 (3), 259–276.



## Stereological formulas

### VERTICAL SPATIAL GRID

<b>Estimated volume</b>	$\hat{V} = a \cdot d_z \cdot \sum_{i=1}^m P_i$	<p><math>a</math> Area associated with point</p> <p><math>d_z</math> Distance between planes</p> <p><math>m</math> Number of scanning planes</p> <p><math>\sum P</math> Intersections with points</p>
<b>Area associated with point</b>	$a = \frac{w^2}{2\pi}$	<p><math>w</math> Horizontal width</p>
<b>Estimated surface area</b>	$\hat{S} = 2 \cdot \frac{a \cdot d_z}{l + \frac{4}{\pi} \cdot d_z} \cdot (I_{xy} + I_{xz})$	<p><math>a</math> Area associated with point</p> <p><math>d_z</math> Distance between planes</p> <p><math>l</math> Length of cycloid</p> <p><math>I_{xy}</math> X,Y intersections</p> <p><math>I_{xz}</math> X,Z intersections</p>
<b>Length of cycloid</b>	$l = \frac{2w}{\pi}$	<p><math>w</math> Horizontal width</p>
<b>X,Y intersections</b>	$I_{xy} = \sum_{i=1}^m I_{xy,i}$	<p><math>m</math> Number of scanning planes</p>
<b>X,Z intersections</b>	$I_{xz} = 2 \cdot \left( \sum_{i=1}^m P_i - \sum_{i=1}^{m-1} P_{i,i+1} \right)$	<p><math>m</math> Number of scanning planes</p> <p><math>\sum P</math> Intersections with points</p>

### References

Cruz-Orive, L. M., Howard, C.V. (1995). Estimation of individual feature surface area with the vertical spatial grid. *Journal of Microscopy*, 178 (2), 146–151.



## Stereological formulas

### WEIBEL

<b>Surface area per unit volume (<math>S_v</math>)</b>	$S_v = 2 \cdot \left( \frac{I}{\frac{l}{2} \cdot P} \right)$	<p><math>I</math> Intersections (triangular markers) <math>P</math> Points (end points circular markers) <math>l</math> Length of each line</p> <p><i>Note: We use <math>l/2</math> for the length represented at each point since there are two end points per line.</i></p>
--	--	---

### References

Weibel, E.R., Kistler, G.S., & Scherle, W.F. (1966). Practical stereological methods for morphometric cytology. *The Journal of cell biology*, 30 (1), 23–38.