

Stereo Investigator Stereological Formulas

Area Fraction Fractionator	2
Cavalieri Estimator	3
Combined Slope Intercept	5
Connectivity Assay	6
Cycloids for L_v	7
Cycloids for S_v	9
Fractionator & Optical Fractionator	12
Isotropic Fakir	15
Isotropic Virtual Planes	16
IUR Planes Optical Fractionator	19
L-Cycloid Optical Fractionator	20
Merz	21
Nucleator	22
Optical Fractionator & Fractionator	23
Optical Rotator	26
Petrimetrics	29
Physical Fractionator	30
Planar Rotator	32
Point Sampled Intercept	33
Size Distribution	34
Spaceballs	36
Surface weighted star volume	38
Surfactor	39
S_v -Cycloid Fractionator	40
Vertical Spatial Grid	41
Weibel	42

Area Fraction Fractionator

Estimated volume fraction	$\hat{V}_v(y, ref) = \frac{\sum_{i=1}^m P(Y)_i}{\sum_{i=1}^m P(ref)_i}$	$P(ref)$ Points hitting reference volume Y Sub-region $P(Y)$ Points hitting sub-region
Estimated area	$\hat{A} = \frac{1}{asf} * a(p) * P(Y_i)$	asf Area sampling fraction $a(p)$ Area associated with a point

References

Howard, C. V., & Reed, M. G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (pp. 170–172). Milton Park, England: BIOS Scientific Publishers.

Cavalieri Estimator

Area associated with a point	$A_p = g^2$	g^2 : Grid area
Volume associated with a point	$V_p = g^2 m \bar{t}$	g^2 : Grid area m : Section evaluation interval $\bar{t}$: Mean section cut thickness
Estimated volume	$\hat{V} = A_p m' \bar{t} \left(\sum_{i=1}^n P_i \right)$	A_p : Area associated with a point g : Grid size m' : Section evaluation interval $\bar{t}$: Mean section cut thickness P_i : Points counted on grid
Estimated volume corrected for over-projection	$[v] = t \cdot \left(k \cdot \sum_{j=1}^g a'_j - \max(a') \right)$	t : Section cut thickness k : Correction factor g : Grid size a' : Projected area
Variance of systematic random sampling	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n P_i^2, \quad B = \sum_{i=1}^{n-1} P_i P_{i+1}, \quad C = \sum_{i=1}^{n-2} P_i P_{i+2}$ n: Number of sections s^2: Variance due to noise $s^2 = 0.0724 \left(\frac{b}{\sqrt{a}} \right) \sqrt{n \sum_{i=1}^n P_i}$ where $\frac{b}{\sqrt{a}}$ is the shape factor m: Smoothness class of sampled function

Cavalieri Estimator

Coefficient of Error	$CE = \frac{\sqrt{TotalVar}}{\sum_{i=1}^n P_i}$	<i>TotalVar</i> : Total variance of the estimated volume $TotalVar = s^2 + VAR_{SRS}$ <i>n</i>: Number of sections <i>P_i</i>: Points counted on grid
-----------------------------	---	---

References

- García-Fiñana, M., Cruz-Orive, L.M., Mackay, C.E., Pakkenberg, B. & Roberts, N. (2003). [Comparison of MR imaging against physical sectioning to estimate the volume of human cerebral compartments.](#) *Neuroimage*, 18 (2), 505–516.
- Gundersen, H. J. G., & Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction.](#) *Journal of Microscopy*, 147 (3), 229–263.
- Howard, C. V., & Reed, M.G. (2005). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (Chapter 3). New York: Garland Science/BIOS Scientific Publishers.

Combined Slope Intercept

Profile area	$a = a(p) \cdot \sum P$	$a(p)$: Area associated with a point $\sum P$: Number of points
Profile boundary	$b = \frac{\pi}{2} d \cdot \sum I$	d : Distance between points $\sum I$: Number of intersections

This method is based on the principles described in the following:

Miles, R.E., & Davy, P. (1976). Precise and general conditions for the validity of a comprehensive set of stereological fundamental formulae. *Journal of Microscopy*, 107(3), 211-226.

Connectivity Assay

Euler number (X_3):	$X_3 = I + H - B$	I : Total island markers H : Total hole markers B : Total bridge markers
Number of Alveoli (N_{alv})	$N_{alv} = -X_3$	X_3 : Euler number
Sum counting frame volumes (V)	$V = h \cdot n \cdot a$	h : Disector height n : Number of disectors a : Area counting frame
Numerical density of alveoli (N_v)	$N_v = \frac{N_{alv}}{V}$	N_{alv} : Number of alveoli V : Sum counting frame volumes

References

Ochs, M., Nyengaard, J.R., Jung, A., Knudsen, L., Voigt, M., Wahlers, T., Richter, J., Gundersen, H.J.G. (2004). [The number of alveoli in the human lung](#). *American journal of respiratory and critical care medicine*, 169 (1), 120–124.

Cycloids for L_v

Area associated with a point	$A_p = g^2$	g^2 : Grid area
Volume associated with a point	$V_p = g^2 m \bar{t}$	m : Section evaluation interval $\bar{t}$: Mean section cut thickness
Length per unit volume	$L_V = 2 \frac{[\bar{I}_L^c]_{prj}}{\Delta}$ $L_V = \frac{2}{\Delta} \cdot \frac{(\bar{I}_c^{cyc})_{prj}}{\bar{p} \cdot \left(\frac{l}{p}\right)} = \frac{2}{\Delta} \left(\frac{p}{l}\right) \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	$[\bar{I}_L^c]_{prj}$: Number of counting frames Δ : Section cut thickness I_i : Intercepts P_i : Test points $(\bar{I}_c^{cyc})_{prj}$: Average number of intersections of projected images p/l : Test points per unit length of cycloid
Estimated volume	$\hat{V} = m \Delta \left(\frac{a}{p}\right) \sum_{i=1}^n P_i$	m : Sampling fractions Δ : Section cut thickness a : Area p : Number of test points
Estimated length	$\hat{L} = 2 \left(\frac{a}{l}\right) m \sum_{i=1}^n I_i$	a : Area l : Line length m : Sampling fractions I_i : Intercepts
Coefficient of Error for line length	$CE(\hat{L} L) = \frac{\sqrt{VAR_{SRS}}}{\sum_{i=1}^n I_i}$	VAR_{SRS} : Variance of systematic random sampling $\hat{L} L$: Estimated length per length I_i : Intercepts

Cycloids for L_v

Coefficient of Error for length density	$CE(L_v) = \sqrt{\frac{n}{n-1} \left(\frac{\sum_{i=1}^n l_i^2}{\sum_{i=1}^n l_i * \sum_{i=1}^n l_i} + \frac{\sum_{i=1}^n P_i^2}{\sum_{i=1}^n P_i * \sum_{i=1}^n P_i} - 2 * \frac{\sum_{i=1}^n l_i P_i}{\sum_{i=1}^n l_i * \sum_{i=1}^n P_i} \right)}$	l_i : Intercepts P_i : Test points n : Number of probes
Variance of systematic random sampling (VAR_{SRS}):	$VAR_{SRS} = \frac{3g_0 - 4g_1 + g_2}{12}$	g : grid size $g_k = \sum_{i=1}^{n-k} L_i L_{i+k}$ L_i : Line length at section i

References

- Artacho-Pérula, E., & Roldán-Villalobos, R. (1995). [Estimation of capillary length density in skeletal muscle by unbiased stereological methods: I. Use of vertical slices of known thickness.](#) *The Anatomical Record*, 241 (3), 337–344.
- Gokhale, A. M. (1990). [Unbiased estimation of curve length in 3-D using vertical slices.](#) *Journal of Microscopy* 159 (2), 133–141.
- Howard, C. V., & Reed, M. G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy*. Milton Park, England: BIOS Scientific Publishers (pp. 170–172).

Cycloids for Sv

Area associated with a point	$A_p = g^2$	g^2 : Grid area
Volume associated with a point	$V_p = g^2 m \bar{t}$	m : Section evaluation interval $\bar{t}$: Mean section cut thickness
Estimated surface area per unit volume	$est S_v = 2 * \left(\frac{2p}{l}\right) * \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	$\frac{p}{l}$: Points per unit length of cycloid I_i : Intercepts with cycloids P_i : Point counts
Estimated volume	$\hat{V} = m \bar{t} * \frac{a}{p} * \sum_{i=1}^m P_i$	m : Section evaluation interval $\bar{t}$: Mean section cut thickness $\frac{a}{p}$: Area associated with each point P_i : Point counts
Estimated surface area	$\hat{S} = 2 * \frac{a}{l} * m \bar{t} * \sum_{i=1}^m I_i$	$\frac{a}{l}$: Area per unit length m : Section evaluation interval $\bar{t}$: Mean section cut thickness I_i : Intercepts with cycloids
Coefficient of Error for estimated surface	$CE(\hat{S} S) = \frac{\sqrt{VAR_{SRS}}}{\sum_{i=1}^n I_i}$	I_i : Intercepts with cycloids
Variance due to systematic random sampling	$VAR_{SRS} = \frac{3g_0 - 4g_1 + g_2}{12}$ $g_k = \sum_{i=1}^{n-k} S_i * S_{i+k}$	g_0, g_1, g_2 : covariogram values S_i : Surface area at section i

Coefficient of Error for surface density	$CE(S_V) = \sqrt{\frac{n}{n-1} \left(\frac{\sum_{i=1}^n I_i^2}{\sum_{i=1}^n I_i \sum_{i=1}^n I_i} + \frac{\sum_{i=1}^n P_i^2}{\sum_{i=1}^n P_i \sum_{i=1}^n P_i} - 2 \frac{\sum_{i=1}^n I_i P_i}{\sum_{i=1}^n I_i \sum_{i=1}^n P_i} \right)}$	n: Number of measurements I_i : Intercepts P_i : Point counts
--	--	---

References

- Baddeley, A. J., Gundersen, H.J.G., & Cruz-Orive, L.M. (1998). [Estimation of surface area from vertical sections.](#) *Journal of Microscopy* 142 (3), 259–276.
- Howard, C. V., & Reed, M.G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (170–172). Milton Park, England: BIOS Scientific Publishers.

Discrete Vertical Rotator

<u>Estimated volume</u>	$\text{est } v = \frac{\pi}{n} * a_p * \sum_{i=1}^n P_i \cdot D_i$	<u>n</u> : Number of centriolar sections <u>a_p</u> : Area associated with each point <u>P_i</u> : Number of points in each class <u>D_i</u> : Distance of class from central axis
--------------------------------	--	---

-

References

Mironov, A. A. (1998). [Estimation of subcellular organelle volume from ultrathin sections through centrioles with a discretized version of the vertical rotator](#). *Journal of microscopy*, 192 (1), 29–36.

Fractionator & Optical Fractionator

Estimate of total number of particles	<i>Optical Fractionator (3D):</i> $N = \sum Q^- * \frac{t}{h} * \frac{1}{asf} * \frac{1}{ssf}$ <i>Fractionator (2D):</i> $N = \sum Q^- * \frac{1}{asf} * \frac{1}{ssf}$	Q^-: Particles counted t: Section mounted thickness h: Counting frame height asf: Area sampling fraction ssf: Section sampling fraction
Variance due to systematic random sampling (Gundersen)	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- * Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- * Q_{i+2}^-$
Variance due to noise (Gundersen)	$s^2 = \sum_{i=1}^n Q^-$	Q^- : Particles counted
Total variance (Gundersen)	$TotalVar = s^2 + VAR_{SRS}$	VAR_{SRS}: Variance due to systematic random sampling s^2: Variance due to noise
Coefficient of Error (Gundersen)	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$: Total variance s^2: Variance due to noise

Fractionator & Optical Fractionator

Number weighted mean section cut thickness	$\overline{t_Q} = \frac{\sum_{i=1}^m t_i Q_i}{\sum_{i=1}^m Q_i}$	m: Number of sections t_i: Section thickness at site i Q_i: Particles counted
Coefficient of error (Sheaffer)	$CE = \frac{\sqrt{\left(\frac{1}{f} - \frac{1}{F}\right) * s^2}}{\bar{Q}}$	f: Number of counting frames F: Total possible sampling sites s^2: Estimated variance $\bar{Q}$: Average particles counted
Average number of particles (Scheaffer)	$\bar{Q} = \frac{\sum_{i=1}^f Q_i}{f}$	Q_i: Particles counted f: Number of counting frames
Estimated variance (Sheaffer)	$s^2 = \frac{\sum_{i=1}^f (Q_i - \bar{Q})^2}{f - 1}$	f: Number of counting frames Q_i: Particles counted $\bar{Q}$: Average particles counted
Estimated variance of estimated cell population (Sheaffer)	$\frac{C_{fp} F^2 s^2}{f}$	C_{fp}: Finite population correction F: Total possible sampling sites f: Number of counting frames
Estimated variance of mean cell count (Sheaffer)	$\frac{C_{fp} s^2}{f}$	C_{fp}: Finite population correction f: Number of counting frames
Estimated mean coefficient of error (Cruz-Orive)	$est\ Mean\ CE\ (est\ N) = \left[\frac{1}{3n} * \sum_{i=1}^n \left(\frac{Q_{1i} - Q_{2i}}{Q_{1i} Q_{2i}} \right)^2 \right]^{1/2}$	Q_{1i}: Counts in sub-sample 1 Q_{2i}: Counts in sub-sample 2 n: Size of sub-sample

Fractionator & Optical Fractionator

Predicted coefficient of error for estimated population (Schmitz-Hof)	$CE_{pred}(n_F) = \sqrt{\frac{Var(Q_r^-)}{R * (\bar{Q}_r^-)^2}}$ $CE_{pred}(n_F) = \frac{1}{\sqrt{\sum_{r=1}^R Q_r^-}} = \frac{1}{\sqrt{\sum_{s=1}^S Q_s^-}}$	R: Number of counting spaces S: Number of sections Q_s^-: Counts in the "s"-th section Q_r^-: Counts in the "r"-th counting space
---	---	---

References

- Geiser, M., Cruz-Orive, L.M., Hof, V.I., & Gehr, P. (1990). [Assessment of particle retention and clearance in the intrapulmonary conducting airways of hamster lungs with the fractionator](#). *Journal of Microscopy*, 160 (1), 75–88.
- Glaser, E. M., & Wilson, P.D. (1998). [The coefficient of error of optical fractionator population size estimates: a computer simulation comparing three estimators](#). *Journal of Microscopy* 192 (2), 163–171.
- Gundersen, H.G.J., Vedel Jensen, E.B., Kieu, K. & Nielsen, J. (1999). [The efficiency of systematic sampling in stereology—reconsidered](#). *Journal of Microscopy*, 193 (3), 199–211.
- Gundersen, H. J. G., & Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction](#). *Journal of Microscopy*, 147 (3), 229–263.
- Howard, V., & Reed, M. (2005). *Unbiased stereology: three-dimensional measurement in microscopy* (chapter 12). New York: Garland Science/Bios Scientific Publishers.
- Scheaffer, R., Mendenhall, W., Lyman Ott, B., & Gerow, K. (1996). *Elementary survey sampling* (chapter 7). Boston: PWS-Kent Publishing.
- Schmitz, C., & Hof, P.R. (2000). [Recommendations for straightforward and rigorous methods of counting neurons based on a computer simulation approach](#). *Journal of Chemical Neuroanatomy*, 20 (1), 93–114.
- West, M. J., Slomianka, L., & Gundersen, H.J.G. (1991). [Unbiased stereological estimation of the total number of neurons in the subdivisions of the rat hippocampus using the optical fractionator](#). *The Anatomical Record*, 231 (4), 482–497.

Isotropic Fakir

Estimated total surface area	$est S = 2 * \frac{1}{n} * \sum_{i=1}^n \frac{v}{l_i} * l_i$	n : Number of line sets, always set to 3 $\frac{v}{l_i}$: Inverse of the probe per unit volume l_i : Intercepts with test lines
------------------------------	--	--

References

Kubínová, L., & Janacek, J. (1998). [Estimating surface area by the isotropic fakir method from thick slices cut in an arbitrary direction.](#) *Journal of Microscopy*, 191 (2), 201–211.

Isotropic Virtual Planes

Length per unit volume	$L_v = \frac{2p(box)}{a(plane)} * \frac{\sum Q}{\sum p(ref)}$	$p(box)$: Number of corners considered $a(plane)$: Exact sampling area $p(ref)$: Number of corners in region $\sum Q$: Total number of transects
Estimated total length	$L_2 = \frac{1}{ssf} * \frac{1}{asf} * \frac{1}{hsf} * \frac{1}{psd} * 2 \sum Q$ or $L_2 = \frac{1}{ssf} * \frac{dx * dy}{a(box)} * \frac{\bar{t}}{h(box)} * d * 2 \sum Q$ $asf = \frac{a(box)}{dx * dy}$ $hsf = \frac{h(box)}{\bar{t}}$ $psd = \frac{E[a(plane)]}{v(box)}$	ssf: Section sampling fraction asf: Area sampling fraction hsf: Height sampling fraction psd: Probe sampling density $\sum Q$: Total number of transects $a(box)$: Area of sampling box $h(box)$: Depth of sampling box d: Sampling plane separation dx, dy: Distances in x,y directions between systematically sampled fields $\bar{t}$: Average section thickness $a(plane)$: Sampling plane area E: Expected value $v(box)$: Volume of sampling box
Total plane area	$A = \sum_{i=1}^l \sum_{j=1}^s A_{i,j}$	l: Number of layouts s: Number of sampling sites $A_{i,j}$: Plane area inside of each sampling box for each layout

Isotropic Virtual Planes

Coefficient of error (Gundersen)	$CE = \frac{\sqrt{\frac{3(A - \text{Poisson noise}) - 4B + C}{12} + \text{Poisson noise}}}{\sum Q}, m = 0$ $CE = \frac{\sqrt{\frac{3(A - \text{Poisson noise}) - 4B + C}{240} + \text{Poisson noise}}}{\sum Q}, m = 1$	$\text{Poisson noise} = \sum Q_i$ $A = \sum Q_i^2$ $B = \sum (Q_i * Q_{i+1})$ $C = \sum (Q_i * Q_{i+2})$ A, B, C: Covariogram values Q_i : Particles counted
Plane areas	$\vec{V} = (A, B, C)$ Planes: $Ax + By + Cz + D_i$ $D_i = D + (d * i)$ $\text{Box: } \left\{ \left((x, y, z) \left \begin{array}{l} x_0 \leq x \leq x_0 + b_x \\ y_0 \leq y \leq y_0 + b_y \\ z_0 \leq z \leq z_0 + b_z \end{array} \right. \right) \right\}$ $\text{Area: } (\text{box} \cap \text{planes})$	A, B, C, D: Given constants d: Distance between planes i: Integer x_0, y_0, z_0: Vertex of a sampling box b_x, b_y, b_z: Dimensions of a sampling box
Average number of counts	$\bar{Q} = \frac{\sum_{i=1}^p \sum_{j=1}^{l_j} Q_{ij}}{\sum_{j=1}^p l_j}$	p: Number of probes l_j: Number of layouts in each probe Q_{ij}: Number of counts in each probe and layout

**Total corners of sampling boxes
inside the region of interest**

$$C = \sum_{i=1}^p C_i$$

p : Number of probes
 C_i : Number of sampling boxes inside the region of interest

References

Larsen, J. O., Gundersen, H. J. G., & Nielsen, J. (1998). [Global spatial sampling with isotropic virtual planes: estimators of length density and total length in thick, arbitrarily orientated sections.](#) *Journal of Microscopy*, 191 (3), 238-248.

IUR Planes Optical Fractionator

Estimated length	$est L = 2 * \frac{a}{l} \sum I * \frac{1}{ssf} * \frac{1}{asf} * \frac{t}{h}$	$\frac{a}{l}$: Area per unit length of test line $\sum I$: Number of intersections <i>ssf</i> : Section sampling fraction <i>asf</i> : Area sampling fraction <i>t</i> : Section cut thickness <i>h</i> : Height of counting frame
Area sampling fraction	$asf = \frac{area(Frame)}{area(x,y step)} = \frac{x_{CF} * Y_{CF}}{x_{step} * y_{step}}$	X_{CF}, Y_{CF} : X,Y dimensions of counting frame X_{step}, Y_{step} : X,Y dimensions of grid <i>area(Frame)</i> : Area of counting frame <i>area(x,y,step)</i> : Area of grid

L-Cycloid Optical Fractionator

Estimated length of lineal structure	$est L = 2 * \frac{a}{l} * \sum I * \frac{1}{ssf} * \frac{1}{asf} * \frac{t}{h}$	$\frac{a}{l}$: Area per unit cycloid length $\sum I$: Total intercepts ssf: Section sampling fraction asf: Area sampling fraction t: Section cut thickness h: Height of counting frame
Area sampling fraction	$asf = \frac{area(Frame)}{area(x,y step)} = \frac{x_{CF} * y_{CF}}{x_{step} * y_{step}}$	$area(Frame)$: Area of counting frame $area(x,y step)$: Area of grid x_{CF}, y_{CF}: x,y dimensions of counting frame x_{step}, y_{step}: x,y dimensions of grid

References

Stocks, E. A, McArthur, ., J. C. , Griffen, J. W. & Mouton, P. R. (1996). [An unbiased method for estimation of total epidermal nerve fiber length.](#) *Journal of Neurocytology*, 25 (1), 637–644.

Length of semi-circle	$L = \frac{1}{2} \pi d$	d : Circle diameter
Surface area per unit volume	$S_v = \frac{2 \sum I}{\frac{l}{p} \sum P}$	I : Number of intercepts l/p : Length of semi-circle per point P : Number of points

References

- Howard, C. V., & Reed, M. G. (2005). *Unbiased stereology*. New York: Garland Science (prev. BIOS Scientific Publishers).
- Fritsch, R. S., Weibel, E.R. (1981). Stereological Methods, Vol. 1: Practical Methods for Biological Morphometry. *Zeitschrift für allgemeine Mikrobiologie*, 21(8), 630–630.

Nucleator

Area estimation	$a = \pi l^2$	l: Mean length of rays
Volume estimation	$V_n = \frac{4\pi}{3} l_n^3$	l: Mean length of rays
Estimated coefficient of error	$estCV(R) = \frac{\sqrt{\frac{1}{n-1} \sum_{i=1}^n (R_i - \bar{R})^2}}{\bar{R}}$	n : Number of nucleator estimates R _i : Area/volume estimate for each sampling site
Average area/volume estimate	$\bar{R} = \frac{1}{n} \sum_{i=1}^n R_i$	n: Number of nucleator estimates R _i : Area/volume estimate for each sampling site
Relative efficiency	$CE_n(R) = \frac{CV(R)}{\sqrt{n}}$	CV(R): Estimated coefficient of variation n: Number of nucleator estimates
Geometric mean of the area/volume estimate	$e^{\frac{1}{n} \sum_{i=1}^n \ln R_i}$	n: Number of nucleator estimates R _i : Area/volume estimate for each sampling site

References

Gundersen, H. J. G. (1988). [The nucleator](#). *Journal of Microscopy*, 151 (1), 3–21.

Optical Fractionator & Fractionator

Estimate of total number of particles	<i>Optical Fractionator (3D):</i> $N = \sum Q^- * \frac{t}{h} * \frac{1}{asf} * \frac{1}{ssf}$ <i>Fractionator (2D):</i> $N = \sum Q^- * \frac{1}{asf} * \frac{1}{ssf}$	Q^-: Particles counted t: Section mounted thickness h: Counting frame height asf: Area sampling fraction ssf: Section sampling fraction
Variance due to systematic random sampling (Gundersen)	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- * Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- * Q_{i+2}^-$
Variance due to noise (Gundersen)	$s^2 = \sum_{i=1}^n Q^-$	Q^- : Particles counted
Total variance (Gundersen)	$TotalVar = s^2 + VAR_{SRS}$	VAR_{SRS}: Variance due to systematic random sampling s^2: Variance due to noise
Coefficient of Error (Gundersen)	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$: Total variance s^2: Variance due to noise

Optical Fractionator & Fractionator

Number weighted mean section cut thickness	$\overline{t_Q} = \frac{\sum_{i=1}^m t_i Q_i}{\sum_{i=1}^m Q_i}$	m : Number of sections t_i : Section thickness at site i Q_i : Particles counted
Coefficient of error (Sheaffer)	$CE = \frac{\sqrt{\left(\frac{1}{f} - \frac{1}{F}\right) * s^2}}{\bar{Q}}$	f : Number of counting frames F : Total possible sampling sites s^2 : Estimated variance $\bar{Q}$: Average particles counted
Average number of particles (Scheaffer)	$\bar{Q} = \frac{\sum_{i=1}^f Q_i}{f}$	Q_i : Particles counted f : Number of counting frames
Estimated variance (Sheaffer)	$s^2 = \frac{\sum_{i=1}^f (Q_i - \bar{Q})^2}{f - 1}$	f : Number of counting frames Q_i : Particles counted $\bar{Q}$: Average particles counted
Estimated variance of estimated cell population (Sheaffer)	$\frac{C_{fp} F^2 s^2}{f}$	C_{fp} : Finite population correction F : Total possible sampling sites f : Number of counting frames
Estimated variance of mean cell count (Sheaffer)	$\frac{C_{fp} s^2}{f}$	C_{fp} : Finite population correction f : Number of counting frames
Estimated mean coefficient of error (Cruz-Orive)	$est\ Mean\ CE\ (est\ N) = \left[\frac{1}{3n} * \sum_{i=1}^n \left(\frac{Q_{1i} - Q_{2i}}{Q_{1i} Q_{2i}} \right)^2 \right]^{1/2}$	Q_{1i} : Counts in sub-sample 1 Q_{2i} : Counts in sub-sample 2 n : Size of sub-sample

Optical Fractionator & Fractionator

Predicted coefficient of error for estimated population (Schmitz-Hof)	$CE_{pred}(n_F) = \sqrt{\frac{Var(Q_r^-)}{R * (Q_r^-)^2}}$ $CE_{pred}(n_F) = \frac{1}{\sqrt{\sum_{r=1}^R Q_r^-}} = \frac{1}{\sqrt{\sum_{s=1}^S Q_s^-}}$	R: Number of counting spaces S: Number of sections Q_s^-: Counts in the "s"-th section Q_r^-: Counts in the "r"-th counting space
---	---	---

References

- Geiser, M., Cruz-Orive, L.M., Hof, V.I., & Gehr, P. (1990). [Assessment of particle retention and clearance in the intrapulmonary conducting airways of hamster lungs with the fractionator](#). *Journal of Microscopy*, 160 (1), 75–88.
- Glaser, E. M., & Wilson, P.D. (1998). [The coefficient of error of optical fractionator population size estimates: a computer simulation comparing three estimators](#). *Journal of Microscopy* 192 (2), 163–171.
- Gundersen, H.G.J., Vedel Jensen, E.B., Kieu, K. & Nielsen, J. (1999). [The efficiency of systematic sampling in stereology—reconsidered](#). *Journal of Microscopy*, 193 (3), 199–211.
- Gundersen, H. J. G., & Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction](#). *Journal of Microscopy*, 147 (3), 229–263.
- Howard, V., & Reed, M. (2005). *Unbiased stereology: three-dimensional measurement in microscopy* (chapter 12). New York: Garland Science/Bios Scientific Publishers.
- Scheaffer, R., Mendenhall, W., Lyman Ott, B., & Gerow, K. (1996). *Elementary survey sampling* (chapter 7). Boston: PWS-Kent Publishing.
- Schmitz, C., & Hof, P.R. (2000). [Recommendations for straightforward and rigorous methods of counting neurons based on a computer simulation approach](#). *Journal of Chemical Neuroanatomy*, 20 (1), 93–114.
- West, M. J., Slomianka, L., & Gundersen, H.J.G. (1991). [Unbiased stereological estimation of the total number of neurons in the subdivisions of the rat hippocampus using the optical fractionator](#). *The Anatomical Record*, 231 (4), 482–497.

Optical Rotator


Volume of particle	$\hat{v} = a \sum_i^{+/-} g(P_i)$ $a = k * h$	a: Reciprocal line density k : Length of slice h : Systematic spacing
For vertical slabs and lines parallel to vertical axis	$g(p) = \left\{ \frac{\frac{d_1}{2} \arcsin\left(\frac{t}{d_2}\right)}{\arcsin\left(\frac{t}{d_2}\right)} \middle \begin{matrix} d_2 < t \\ t \leq d_2 \end{matrix} \right\}$	d1: Distance along test line d2: Distance from origin to test line t: Thickness of optical slice z: Distance in z from intercept to origin
For vertical slabs and lines perpendicular to vertical axis	$g(P) = \left\{ \frac{\frac{d_1}{2} \arcsin\left(\frac{t}{\sqrt{d_1^2 + z^2}}\right)}{f(0)} \middle \begin{matrix} \sqrt{d_1^2 + z^2} < t \\ z < t \leq \sqrt{d_1^2 + z^2} \\ t \leq z \end{matrix} \right\}$ $f(x) = x + \frac{\pi}{2} \int_x^{d_1} \frac{1}{\arcsin\left(\frac{t}{\sqrt{u^2 + z^2}}\right)} * du$	d1: Distance along test line t: Thickness of optical slice z: Distance in z from intercept to origin

For isotropic slabs	$g(P) = \left\{ \begin{array}{l} \frac{1}{2t} [h(t, d_2) + k(t, d_1, d_2, d_3)] \\ \frac{1}{2t} \left[d_1 d_3 + d_2^2 \log \left(\frac{d_1 + d_3}{t + \sqrt{t^2 - d_2^2}} \right) \right] \end{array} \right\} \begin{array}{l} d_3 < t \\ d_2 < t \leq d_3 \\ t \leq d_2 \end{array}$ $h(t, d) = t^2 \sqrt{1 - \frac{d^2}{t^2}}$ $k(t, d_1, d_2, d_3) = d_1 d_3 + d_2^2 \log \left(\frac{d_1 + d_3}{t + \sqrt{t^2 - d_2^2}} \right)$	 <i>d1</i>: Distance along test line <i>d2</i>: Distance from origin to test line <i>d3</i>: Distance from intercept to origin <i>t</i>: Thickness of optical slice
Estimated surface area	$\hat{s} = a \sum_j l_j g(l_j)$ $g(l) = \left\{ \pi * \frac{2}{\arcsin(\frac{t}{d_2})} \begin{array}{l} 1 \\ t \leq d_2 \end{array} \right\}$	<i>a</i>: Reciprocal line density <i>l_j</i>: Number of intersections between grid line and cell boundary <i>d₂</i>: Distance from origin to test line <i>t</i>: Thickness of optical slice

References

Tandrup, T., Gundersen, H.J.G., & Vedel Jensen, E.B. (1997). [The optical rotator](#). *Journal of microscopy* 186 (2), 108–120.

Petrimetrics

Total length	$\hat{L} = \frac{\pi}{2} * \frac{a}{l} * \frac{1}{asf} * \sum I$ $\hat{L} = d * \frac{1}{asf} * \sum I$	$a/l = 2d/\pi$: Grid constant ($2d/\pi$ units or ratio of area to length of semi-circle probe) asf: Area fraction (ratio of area of counting frame to grid-step) I: Number of intersections counted $d = 2 * \text{Merz-radius}$: The Merz-radius refers to the radius of the semi-circle used to probe. 
----------------------------	--	---

References

Howard, C. V., & Reed, M. G. (2005). *Unbiased stereology*. New York: Garland Science (prev. BIOS Scientific Publishers).

Physical Fractionator

Total number of particles	$N = \sum Q^- * \frac{1}{asf} * \frac{1}{ssf}$	Q^- : Particles counted asf : Area sampling fraction ssf : Section sampling fraction
Variance due to noise	$s^2 = \sum_{i=1}^n Q^-$	n : Number of sections used Q^- : Particles counted
Variance due to systematic random sampling	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$ $A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- * Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- * Q_{i+2}^-$	s^2 : Variance due to noise m : Smoothness class of sampled function
Total variance	$TotalVar = s^2 + VAR_{SRS}$	s^2 : Variance due to noise VAR_{SRS} : Variance due to systematic random sampling
Coefficient of error	$CE = \frac{\sqrt{TotalVar}}{s^2}$	$TotalVar$: Total variance s^2 : Variance due to noise

References

Gundersen, H.J.G. (1986). [Stereology of arbitrary particles.](#) *Journal of Microscopy*, 143 (1), 3–45.

Sterio, D. C. (1984). [The unbiased estimation of number and sizes of arbitrary particles using the disector.](#) *Journal of Microscopy* 134 (2), 127–136.

Planar Rotator

Volume for the isotropic planar rotator	$V = 2t \sum_i g_i$	t : Separation between test lines g_i : Isotropic planar rotator function
Volume for the vertical planar rotator	$V = \pi t \sum_i l_i^2$	t : Separation between test lines l_i : Intercept length along a test line
Isotropic planar rotator function	$g_i(l) = l \sqrt{l^2 + a_i^2} + a_i^2 \ln \left[\frac{l}{a_i} + \sqrt{\left(\frac{l}{a_i}\right)^2 + 1} \right]$ $g_{i+} = \sum_{j \text{ even}} g_i(l_{ij+}) - \sum_{j \text{ odd}} g_i(l_{ij+})$ $g_{i-} = \sum_{j \text{ even}} g_i(l_{ij-}) - \sum_{j \text{ odd}} g_i(l_{ij-})$ $g_i = \frac{1}{2}(g_{i+} + g_{i-})$	l : Intercept length along a test line a_i : Distance from origin to test line j : Number of grid lines l_{ij} : Number of intersections between the j -th grid line and the cell boundary
	$l_{i+}^2 = \sum_{j \text{ even}} l_{ij+}^2 - \sum_{j \text{ odd}} l_{ij+}^2$ $l_{i-}^2 = \sum_{j \text{ even}} l_{ij-}^2 - \sum_{j \text{ odd}} l_{ij-}^2$ $l_i^2 = \frac{1}{2}(l_{i+}^2 + l_{i-}^2)$	l : Intercept length along a test line l_{ij} : Number of intersections between the j -th grid line and the cell boundary

References

Jensen, E.B., Gundersen. H.J.G. (1993). [The rotator](#). *Journal of Microscopy* 170 (1), 35–44.

Point Sampled Intercept

Volume based on intercept length	$\hat{V}_V = \frac{\pi}{3} \bar{l}_0^3 = \frac{\pi}{3n} \sum_{i=1}^n l_{0,i}^3$	n : Number of intercepts l : Intercept length
Volume weighted mean volume	$\bar{v}_V = \frac{\sum_{i=1}^n \bar{l}_0^3}{n} * \frac{\pi}{3}$	n : Number of intercepts l : Intercept length
Coefficient of error	$CE(\bar{l}_0^3) = \sqrt{\frac{\sum_{i=1}^n (\bar{l}_0^3)^2}{(\sum_{i=1}^n \bar{l}_0^3)^2} - \frac{1}{n}}$	n : Number of intercepts l : Intercept length
Coefficient of variance	$CV(\bar{l}_0^3) = CE(\bar{v}_V) * \sqrt{n}$	l : Intercept length $\bar{v}_V$: Volume weighted mean volume n : Number of intercepts
Variance	$Variance V(v) = \left[\frac{\pi}{3} * SD(\bar{l}_0^3) \right] = [CV(\bar{l}_0^3) * \bar{v}_V]$	l : Intercept length $CV(\bar{l}_0^3)$: Coefficient of variance $\bar{v}_V$: Volume weighted mean volume

References

- Gundersen, H.J.G., Jensen. E.B. (1985). [Stereological Estimation of the Volume-Weighted Mean Volume of Arbitrary Particles Observed on Random Sections](#). *Journal of Microscopy*, 138, 127–142.
- Sørensen, F.B. (1991). [Stereological estimation of the mean and variance of nuclear volume from vertical sections](#). *Journal of Microscopy*, 162, (2), 203–229.

Size Distribution

Volume-weighted mean particle volume	$\bar{v}_V = \bar{v}_N * [1 + CV_N^2(v)]$	$\bar{v}_N$: Number-weighted mean volume $CV_N(v)$: Coefficient of variation
Number-weighted mean volume	$\bar{v}_N = \frac{\sum S}{\sum R}$ $S = Q * R$ $T = Q^2 * R$	R : Number of contours Q : Number of points per contour
Variance	$Var_N(v) = \frac{\left[\sum T - \frac{(\sum S)^2}{\sum R} \right] * v(p)^2}{\sum R - 1}$	$v(p)$: Volume associated with a point
Standard deviation	$SD_N(v) = \sqrt{Var_N(v)}$ $SD_N(v) = \sqrt{\bar{v}_N * (\bar{v}_V - \bar{v}_N)}$	$\bar{v}_N$: Number-weighted mean volume $\bar{v}_V$: Volume-weighted particle volume
Coefficient of variation	$CV_N(v) = \frac{SD_N(v)}{\bar{v}_N}$ $CV_N(v) = \sqrt{\frac{\bar{v}_V - \bar{v}_N}{\bar{v}_N}}$	SD_N : Standard deviation $\bar{v}_N$: Number-weighted mean volume $\bar{v}_V$: Volume-weighted particle volume
Coefficient of error	$CE_N(v) = \frac{CV_N(v)}{\sqrt{R}}$	$CV_N(v)$: Coefficient of variation R : Number of contours

Size Distribution

References

Sørensen, F.B. (1991). [Stereological estimation of the mean and variance of nuclear volume from vertical sections](#). *Journal of Microscopy*, 162 (2), 203–229.

Spaceballs

Length estimate	$L = 2 * \left(\sum_{i=1}^n Q_i \right) * \frac{v}{a} * \frac{1}{ssf}$ i This equation is equivalent to the equation used by Mouton et al.	n: Number of sections used Q_i: Intersections counted v: Volume (grid X * grid Y * section thickness) a: Surface area of the sphere ssf: Section sampling fraction
Variance due to noise	$s^2 = \sum_{i=1}^n Q_i$	Q_i: Intersections counted
Variance due to systematic random sampling	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$ $A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- * Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- * Q_{i+2}^-$	s^2: Variance due to noise m: Smoothness class of sampled function
Total variance	$TotalVar = s^2 + VAR_{SRS}$	s^2: Variance due to noise VAR_{SRS}: Variance due to systematic random sampling

Coefficient of error	$CE = \frac{\sqrt{TotalVar}}{s^2}$	<i>TotalVar</i> : Total variance s^2 : Variance due to noise
-----------------------------	------------------------------------	---

References

- Gundersen, H.J.G., Vedel Jensen, E.B., Kieu, K. & Nielsen, J. (1999). [The efficiency of systematic sampling in stereology—reconsidered](#). *Journal of Microscopy*, 193 (3) 199–211.
- Gundersen, H. J. G., Vedel Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction](#). *Journal of Microscopy*, 147 (3), 229–263.
- Mouton, P. R., Gokhale, A.M., Ward, N.L., & West, M.J. (2002). [Stereological length estimation using spherical probes](#). *Journal of Microscopy*, 206 (1), 54–64.

Surface weighted star volume

Surface-weighted star volume	$\hat{v}_s^* = \frac{2\pi}{3} * \bar{l}_1^3$ $\hat{v}_s^* = \frac{2\pi}{3} * \frac{\sum_{i=1}^n \sum_{j=1}^{m_i} l_{1,(i,j)}^3}{\sum_{i=1}^n m_i}$	l : Intercept length m_i : Number of intercepts n : Number of probes
Sum of cubed intercepts in probe	$y_i = \sum_{j=1}^{m_i} l_{1,(i,j)}^3$	l : Intercept length m_i : Number of intercepts
Coefficient of error	$CE[\hat{v}_s^*] = \left[\frac{n}{n-1} \left\{ \frac{\sum y_i^2}{\sum y_i \sum y_i} + \frac{\sum m_i^2}{\sum m_i \sum m_i} - 2 \frac{\sum m_i * y_i}{\sum m_i \sum y_i} \right\} \right]^{1/2}$	n : Number of probes y_i : Sum of cubed intercepts in probe m_i : Number of intercepts

References

Reed, M. G., Howard, C.V. (1998). [Surface-weighted star volume: concept and estimation](#). *Journal of Microscopy*, 190 (3), 350–356.

Surfactor

Surface area for single-ray designs	$\hat{S} = 4\pi l_0^2 + c(\beta)$	l : Length of intercept β : Angle between test line and surface $c()$: Function of the planar angle
Surface area for multi-ray designs	$\hat{S} = 2\pi \sum_{j=1}^{2r} l_j^2 * c(\beta)$	l : Length of intercept β : Angle between test line and surface $c()$: Function of the planar angle r : Number of test lines
Function of the planar angle	$c(\beta) = 1 + \left[\frac{1}{2} \cos \beta \right] * \left[\frac{\pi}{2} - \sin^{-1} \frac{1 - \cos^2 \beta}{1 + \cos^2 \beta} \right]$	β : Angle between test line and surface

References

Jensen, E.B., Gundersen, H.J.G. (1987). Stereological estimation of surface area of arbitrary particles. *Acta Stereologica*, 6 (3).

Sv-Cycloid Fractionator

Estimated surface area per unit volume	$S_V = 2 \left(\frac{p}{l}\right) * \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	p/l: Ratio of test points to curve length $\sum I_i$: Total intercept points on curve $\sum P_i$: Total test points n: Number of micrographs
--	--	--

References

Baddeley, A. J., Gundersen, H.J.G., & Cruz-Orive, L.M. (1986). [Estimation of surface area from vertical sections.](#) *Journal of Microscopy*, 142 (3), 259–276.

Vertical Spatial Grid

Estimated volume	$\hat{V} = a * d_z * \sum_{i=1}^m P_i$	a : Area associated with point d_z : Distance between Z planes m : Number of scanning planes $\sum P_i$: Intersections with points
Area associated with point	$a = \frac{w^2}{2\pi}$	w : Horizontal width
Estimated surface area	$\hat{S} = 2 * \frac{a * d_z}{l + \frac{4}{\pi} * d_z} * (I_{xy} + I_{xz})$	a : Area associated with point d_z : Distance between Z planes l : Length of cycloid I_{xy} : X,Y intersections I_{xz} : X,Z intersections
Length of cycloid	$l = \frac{2w}{\pi}$	w : Horizontal width
X,Y intersections	$I_{xy} = \sum_{i=1}^m I_{xy,i}$	m : Number of scanning planes
X,Z intersections	$I_{xz} = 2 * \left(\sum_{i=1}^m P_i - \sum_{i=1}^{m-1} P_{i,i+1} \right)$	m : Number of scanning planes $\sum P_i$: Intersections with points

References

Cruz-Orive, L. M., Howard, C.V. (1995). [Estimation of individual feature surface area with the vertical spatial grid](#). *Journal of Microscopy*, 178 (2), 146–151.

Weibel

Surface area per unit volume	$S_V = 2N_i = \frac{2N}{d}$	d : Length of test line N_i : Number of intercepts per total length N : Number of intercepts
-------------------------------------	-----------------------------	--

References

Weibel, E.R., Kistler, G.S., & Scherle, W.F. (1966). [Practical stereological methods for morphometric cytology](#). *The Journal of cell biology*, 30 (1), 23–38.