



Stereological formulas

AREA FRACTION FRACTIONATOR	2
CAVALIERI ESTIMATOR.....	3
CYCLOIDS FOR LV	5
IMAGE VOLUME FRACTIONATOR.....	7
IMAGE VOLUME SPACEBALLS	11
ISOTROPIC FAKIR.....	13
NUCLEATOR	14
OPTICAL ROTATOR.....	15
PLANAR ROTATOR.....	17
SURFACTOR	19

AREA FRACTION FRACTIONATOR

<p>Estimated volume fraction (\hat{V}_v)</p>	$\hat{V}_v(Y, ref) = \frac{\sum_{i=1}^m P(Y)_i}{\sum_{i=1}^m P(ref)_i}$	<p>$P(ref)$ Points hitting reference volume Y Sub-region $P(Y)$ Points hitting sub-region</p>
<p>Estimated area (\hat{A})</p>	$\hat{A} = \frac{1}{asf} \cdot a(p) \cdot P(Y_i)$	<p>asf Area sampling fraction $a(p)$ Area associated with a point</p>

References

Howard, C. V., & Reed, M. G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (pp. 170–172). Milton Park, England: BIOS Scientific Publishers.

CAVALIERI ESTIMATOR

Area associated with a point (A_p)	$A_p = g^2$	g^2 Grid area
Volume associated with a point (V_p)	$V_p = g^2 m \bar{t}$	m Section evaluation interval \bar{t} Mean section cut thickness
Estimated volume (\hat{V})	$\hat{V} = A_p m' \bar{t} \left(\sum_{i=1}^n P_i \right)$	A_p Area associated with a point m' Section evaluation interval \bar{t} Mean section cut thickness P_i Points counted on grid
Estimated volume corrected for over-projection ($[v]$)	$[v] = t \cdot \left(k \cdot \sum_{j=1}^g a'_j - \max(a') \right)$	t Section cut thickness k Correction factor g Grid size a' Projected area
Coefficient of error (CE)	$CE = \frac{\sqrt{TotalVar}}{\sum_{i=1}^n P_i}$	$TotalVar$ Total variance of the estimated volume n Number of sections P_i Points counted on grid $TotalVar = s^2 + VAR_{SRS}$

Cavalieri Estimator (2)

<p>Variance of systematic random sampling (VAR_{SRS})</p>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	<p>m Smoothness class of sampled function s^2 Variance due to noise $A = \sum_{i=1}^n P_i^2$, $B = \sum_{i=1}^{n-1} P_i P_{i+1}$, $C = \sum_{i=1}^{n-2} P_i P_{i+2}$</p> <p>With:</p> <p>$n$: number of sections</p> $s^2 = 0.0724 \left(\frac{b}{\sqrt{a}}\right) \sqrt{n \sum_{i=1}^n P_i}$ <p>$\frac{b}{\sqrt{a}}$ Shape factor</p>
---	--	--

References

García-Fiñana, M., Cruz-Orive, L.M., Mackay, C.E., Pakkenberg, B. & Roberts, N. (2003). [Comparison of MR imaging against physical sectioning to estimate the volume of human cerebral compartments](#). *Neuroimage*, 18 (2), 505–516.

Gundersen, H. J. G., & Jensen, E.B. (1987). [The efficiency of systematic sampling in stereology and its prediction](#). *Journal of Microscopy*, 147 (3), 229–263.

Howard, C. V., & Reed, M.G. (2005). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (Chapter 3). New York: Garland Science/BIOS Scientific Publishers.

CYCLOIDS FOR LV

<p>Area associated with a point (A_p)</p>	$A_p = g^2$	<p>g^2 Grid area</p>
<p>Volume associated with a point (V_p)</p>	$V_p = g^2 m \bar{t}$	<p>g^2 Grid area m Section evaluation interval \bar{t} Mean section cut thickness</p>
<p>Length per unit volume (L_V)</p>	$L_V = 2 \frac{[\bar{I}_L^C]_{prj}}{\Delta}$ $L_V = \frac{2}{\Delta} \cdot \frac{(\bar{I}_c^{cyc})_{prj}}{\bar{p} \cdot \left(\frac{l}{p}\right)} = \frac{2}{\Delta} \left(\frac{p}{l}\right) \frac{\sum_{i=1}^n I_i}{\sum_{i=1}^n P_i}$	<p>$[\bar{I}_L^C]_{prj}$ Number of counting frames Δ Section cut thickness I_i Intercepts P_i Test points $[\bar{I}_C^{cyc}]_{prj}$ Average number of intersections of projected images $\frac{p}{l}$ Test points per unit length of cycloid</p>
<p>Estimated volume (\hat{V})</p>	$\hat{V} = m \Delta \left(\frac{a}{p}\right) \sum_{i=1}^n P_i$	<p>m Sampling fractions Δ Section cut thickness a Area p Number of test points P_i Test points</p>
<p>Estimated length (\hat{L})</p>	$\hat{L} = 2 \left(\frac{a}{l}\right) m \sum_{i=1}^n I_i$	<p>a Area l Line length m Sampling fractions I_i Intercepts</p>

Cycloids for L_v (2)

Coefficient of error for line length	$CE(\hat{L} L) = \frac{\sqrt{VAR_{SRS}}}{\sum_{i=1}^n I_i}$	VAR_{SRS} Variance of systematic random sampling $\hat{L} L$ Estimated length per length I_i Intercepts
Variance of systematic random sampling (VAR_{SRS})	$VAR_{SRS} = \frac{3g_0 - 4g_1 + g_2}{12}$ $g_k = \sum_{i=1}^{n-k} L_i L_{i+k}$	g Grid size L_i Line length at section i
Coefficient of error for length density	$CE(L_v) = \sqrt{\frac{n}{n-1} \left(\frac{\sum_{i=1}^n I_i^2}{\sum_{i=1}^n I_i \sum_{i=1}^n I_i} + \frac{\sum_{i=1}^n P_i^2}{\sum_{i=1}^n P_i \sum_{i=1}^n P_i} - 2 \frac{\sum_{i=1}^n I_i P_i}{\sum_{i=1}^n I_i \sum_{i=1}^n P_i} \right)}$	I_i Intercepts P_i Test points n Number of probes

References

Artacho-Pérula, E., Roldán-Villalobos, R. (1995). Estimation of capillary length density in skeletal muscle by unbiased stereological methods: I. Use of vertical slices of known thickness *The Anatomical Record*, 241 (3), 337-344.

Gokhale, A. M. (1990). Unbiased estimation of curve length in 3-D using vertical slices. *Journal of Microscopy*, 159 (2), 133–141.

Howard, C. V., Reed, M.G. (1998). *Unbiased Stereology, Three-Dimensional Measurement in Microscopy* (pp. 170–172). BIOS Scientific Publishers.

IMAGE VOLUME FRACTIONATOR

<p>Estimate of total number of particles (N)</p>	$N = \sum Q^- \cdot \frac{1}{asf} \cdot \frac{1}{zsf}$	<p>Q^- Particles counted asf Area sampling fraction (counting frame/grid size) zsf Section sampling fraction (disector height/virtual section thickness)</p>
<p>Variance due to systematic random sampling – Gundersen (VAR_{SRS})</p>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m$ $= 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m$ $= 1$	<p>$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-$ s^2 Variance due to noise</p>
<p>Variance due to noise - Gundersen (s^2)</p>	$s^2 = \sum_{i=1}^n Q^-$	<p>Q^- Particles counted n Number of sections used</p>
<p>Total variance – Gundersen ($TotalVar$)</p>	$TotalVar = s^2 + VAR_{SRS}$	<p>VAR_{SRS} Variance due to SRS s^2 Variance due to noise</p>
<p>Coefficient of error – Gundersen (CE)</p>	$CE = \frac{\sqrt{TotalVar}}{s^2}$	<p>$TotalVar$ Total variance s^2 Variance due to noise</p>

Stereological formulas

Image Volume Fractionator (2)

Coefficient of error – Scheaffer (CE)	$CE = \frac{\sqrt{s^2 \left(\frac{1}{f} - \frac{1}{F} \right)}}{\bar{Q}}$	<i>f</i> Number of counting frames <i>F</i> Total possible sampling sites <i>s</i> ² Estimated variance \bar{Q} Average particles counted
Average number of particles – Scheaffer (\bar{Q})	$\bar{Q} = \frac{\sum_{i=1}^f Q_i}{f}$	Q_i Particles counted <i>f</i> Number of counting frames
Estimated variance - Scheaffer (<i>s</i>²)	$s^2 = \frac{\sum_{i=1}^f (Q_i - \bar{Q})^2}{f - 1}$	<i>f</i> Number of counting frames Q_i Particles counted \bar{Q} Average particles counted
Estimated variance of estimated cell population - Scheaffer	$\frac{C_{fp} F^2 s^2}{f}$	C_{fp} Finite population correction s^2 Estimated variance <i>f</i> Number of counting frames <i>F</i> Total possible sampling sites
Estimated variance of mean cell count - Scheaffer	$\frac{C_{fp} s^2}{f}$	C_{fp} Finite population correction s^2 Estimated variance <i>f</i> Number of counting frames

Image Volume Fractionator (3)

Estimated mean coefficient of error – Cruz-Orive (<i>est Mean CE</i>)	$\text{est Mean CE (est } N) = \left[\frac{1}{3n} \cdot \sum_{i=1}^n \left(\frac{Q_{1i} - Q_{2i}}{Q_{1i} + Q_{2i}} \right)^2 \right]^{1/2}$	Q_{1i} Counts in sub-sample 1 Q_{2i} Counts in sub-sample 2 n Size of sub-sample
Predicted coefficient of error for estimated population – Schmitz-Hof (CE_{pred})	$CE_{pred}(n_F) = \sqrt{\frac{\text{Var}(Q_r^-)}{R \cdot (Q_r^-)^2}}$ $CE_{pred}(n_F) = \frac{1}{\sqrt{\sum_{r=1}^R Q_r^-}} = \frac{1}{\sqrt{\sum_{s=1}^S Q_s^-}}$	R Number of counting spaces S Number of sections Q_r^- Counts in the "r"-th counting space Q_s^- Counts in the "s"-th section

References

- Geiser, M., Cruz-Orive, L.M., Hof, V.I., & Gehr, P. (1990) Assessment of particle retention and clearance in the intrapulmonary conducting airways of hamster lungs with the fractionator. *Journal of Microscopy*, 160 (1), 75–88.
- Glaser, E. M., Wilson, P.D. (1998). The coefficient of error of optical fractionator population size estimates: a computer simulation comparing three estimators. *Journal of Microscopy*, 192 (2), 163–171.
- Gundersen, H.J.G., Vedel Jensen, E.B., Kieu, K., & Nielsen, J. (1999). The efficiency of systematic sampling in stereology—reconsidered. *Journal of Microscopy*, 193 (3), 199–211.
- Gundersen, H. J. G., Jensen, E.B. (1987). The efficiency of systematic sampling in stereology and its prediction. *Journal of Microscopy*, 147 (3), 229–263.
- Howard, V., Reed, M. (2005). *Unbiased stereology: three-dimensional measurement in microscopy* (vol. 4, chapter 12). Garland Science/Bios Scientific Publishers.



Stereological formulas

Image Volume Fractionator (4)

Scheaffer, R.L., Ott, L., & Mendenhall, W. (1996). *Elementary survey sampling* (chapter 7). Boston: PWS-Kent.

Schmitz, C., Hof, P.R. (2000). Recommendations for straightforward and rigorous methods of counting neurons based on a computer simulation approach. *Journal of Chemical Neuroanatomy*, 20 (1), 93–114.

West, M. J., Slomianka, L., & Gundersen, H.J.G. (1991). Unbiased stereological estimation of the total number of neurons in the subdivisions of the rat hippocampus using the optical fractionator. *The Anatomical Record*, 231 (4), 482–497.

IMAGE VOLUME SPACEBALLS

<p>Length estimate</p>	$L = 2 \cdot \left(\sum_{i=1}^n Q_i \right) \cdot \frac{v}{a}$ <p><i>This equation does not include the terms F2 (area-fraction) and F3 (thickness-fraction) used by Mouton et al. (equation 2, 2002), but includes that information in v (volume sampled).</i></p>	<p>n Number of sections used Q_i Intersection counted v Volume (grid X * grid Y * section thickness) a Surface area of the sphere</p>
<p>Variance due to noise</p>	$s^2 = \sum_{i=1}^n Q_i$	<p>Q_i Intersection counted</p>
<p>Variance due to systematic random sampling</p>	$VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{12}, m = 0$ $VAR_{SRS} = \frac{3(A - s^2) - 4B + C}{240}, m = 1$	<p>$A = \sum_{i=1}^n (Q_i^-)^2$ $B = \sum_{i=1}^{n-1} Q_i^- Q_{i+1}^-$ $C = \sum_{i=1}^{n-2} Q_i^- Q_{i+2}^-$</p> <p>$s^2$ Variance due to noise m Smoothness class of sampled function</p>
<p>Total variance</p>	$TotalVar = s^2 + VAR_{SRS}$	<p>VAR_{SRS} Variance due to SRS s^2 Variance due to noise</p>



Stereological formulas

Image Volume Spaceballs (2)

Coefficient of error	$CE = \frac{\sqrt{TotalVar}}{s^2}$	<i>TotalVar</i> Total variance <i>s</i> ² Variance due to noise
-----------------------------	------------------------------------	---

References

Mouton, P. R., Gokhale, A.M., Ward, N.L., & West, M.J. (2002). Stereological length estimation using spherical probes. *Journal of Microscopy*, 206 (1), 54–64.



ISOTROPIC FAKIR

Estimated total surface area	$estS = 2 \frac{1}{n} \cdot \sum_{i=1}^n \frac{v}{l_i} \cdot I_i$	n Number of line sets (always set to 3) $\frac{v}{l_i}$ Inverse of the probe per unit volume I_i Intercepts with test lines
-------------------------------------	---	---

References

Kubínová, L., Janacek, J. (1998). Estimating surface area by the isotropic fakir method from thick slices cut in an arbitrary direction. *Journal of Microscopy*, 191 (2), 201–211.

NUCLEATOR

Area estimate	$a = \pi \bar{l}^2$	l Length of rays
Volume estimate	$\bar{v}_N = \frac{4\pi}{3} \bar{l}_n^3$	l Length of rays
Estimated coefficient of error	$\text{est } CV(R) = \frac{\sqrt{\frac{1}{n-1} \sum_{i=1}^n (R_i - \bar{R})^2}}{\bar{R}}$	n Number of nucleator estimates R_i Area/volume estimate for each sampling site
Average area/volume estimate	$\bar{R} = \frac{1}{n} \sum_{i=1}^n R_i$	n Number of nucleator estimates R_i Area/volume estimate for each sampling site
Relative efficiency	$CE_n(R) = \frac{CV(R)}{\sqrt{n}}$	n Number of nucleator estimates $CV(R)$ Estimated coefficient of variation
Geometric mean of area/volume estimate	$e^{\left(\frac{1}{n} \sum_{i=1}^n \ln R_i\right)}$	n Number of nucleator estimates R_i Area/volume estimate for each sampling site

References

Gundersen, H.J.G. (1988). The nucleator. *Journal of Microscopy*, 151 (1), 3–21.

OPTICAL ROTATOR

<p>Volume of particle</p>	$\hat{v} = a \sum_i^{+/-} g(P_i)$	<p>a Reciprocal line density</p> <p>$a = k \cdot h$</p> <p>k Length of slice</p> <p>h Systematic spacing</p>
<p>For vertical slabs and lines parallel to vertical axis</p>	$g(P) = d_1, \text{ if } d_2 < t$ $g(P) = \frac{\frac{\pi}{2} d_1}{\arcsin\left(\frac{t}{d_2}\right)}, \text{ if } t \leq d_2$	<p>d_1 Distance along test line</p> <p>d_2 Distance from origin to test line</p> <p>t $\frac{1}{2}$ thickness of optical slice</p>
<p>For vertical slabs and lines perpendicular to vertical axis</p>	$g(P) = d_1, \text{ if } \sqrt{d_1^2 + z^2} < t$ $g(P) = f\left(\sqrt{t^2 - z^2}\right), \text{ if } z < t \leq \sqrt{d_1^2 + z^2}$ $g(P) = f(0), \text{ if } t \leq z $ $f(x) = x + \frac{\pi}{2} \int_x^{d_1} \frac{1}{\arcsin\left(\frac{t}{\sqrt{u^2 + z^2}}\right)} du$	<p>d_1 Distance along test line</p> <p>t $\frac{1}{2}$ thickness of optical slice</p> <p>z Distance in z from intercept to origin</p>

Stereological formulas

Optical Rotator (2)

<p>For isotropic slabs</p>		$g(P) = d_1, \quad \text{if } d_3 < t$ $g(P) = \frac{1}{2t} [h(t, d_2) + k(t, d_1, d_2, d_3)], \text{ if } d_2 < t \leq d_3$ $h(t, d) = t^2 \sqrt{1 - \frac{d^2}{t^2}}$ $k(t, d_1, d_2, d_3) = d_1 d_3 + d_2^2 \log \left(\frac{d_1 + d_3}{t + \sqrt{t^2 - d_2^2}} \right)$	<p>d_1 Distance along test line d_2 Distance from origin to test line d_3 Distance from intercept to origin t ½ thickness of optical slice</p>
<p>Estimated surface area</p>		$\hat{S} = a \sum_j l_j g(l_j)$ $g(l) = 2, \quad \text{if } d_2 < t$ $g(l) = \pi \cdot \frac{1}{\arcsin\left(\frac{t}{d_2}\right)}, \quad \text{if } t \leq d_2$	<p>a Reciprocal line density l_j Number of intersections between grid line and cell boundary d_2 Distance from origin to test line t ½ thickness of optical slice</p>

References

Tandrup, T., Gundersen, H.J.G., & Vedel Jensen, E.B. (1997). The optical rotator *Journal of microscopy*, 186 (2), 108–120.

PLANAR ROTATOR

<p>Volume for isotropic planar rotator</p>	$V = 2t \sum_i g_i$	<p>t Separation between test lines g_i Isotropic planar rotator function</p>
<p>Volume for vertical planar rotator</p>	$V = \pi t \sum_i l_i^2$	<p>t Separation between test lines l_i Intercept length along a test line</p>
<p>Isotropic planar rotator function</p>	$g_i(l) = l \sqrt{l^2 + a_i^2} + a_i^2 \ln \left[\frac{l}{a_i} + \sqrt{\left(\frac{l}{a_i}\right)^2 + 1} \right]$ $g_{i+} = \sum_{j \text{ even}} g_i(l_{i j+}) - \sum_{j \text{ odd}} g_i(l_{i j+})$ $g_{i-} = \sum_{j \text{ even}} g_i(l_{i j-}) - \sum_{j \text{ odd}} g_i(l_{i j-})$ $g_i = \frac{1}{2} (g_{i+} + g_{i-})$	<p>l Intercept length along a test line a_i Distance from origin to test line j Number of grid lines l_{ij} Number of intersections between the j-th grid line and the cell boundary</p>

Planar Rotator (2)

Isotropic planar rotator function (cont'd)	$l_{i+}^2 = \sum_{j \text{ even}} l_{ij+}^2 - \sum_{j \text{ odd}} l_{ij+}^2$ $l_{i-}^2 = \sum_{j \text{ even}} l_{ij-}^2 - \sum_{j \text{ odd}} l_{ij-}^2$ $l_i^2 = \frac{1}{2}(l_{i+}^2 + l_{i-}^2)$	<p>l Intercept length along a test line</p> <p>a_i Distance from origin to test line</p> <p>j Number of grid lines</p> <p>l_{ij} Number of intersections between the j-th grid line and the cell boundary</p>
---	--	--

References

Jensen Vedel, E.B., Gundersen, H.J.G. (1993). The rotator *Journal of Microscopy*, 170 (1), 35–44.

SURFACTOR

Surface area for single-ray designs	$\hat{S} = 4\pi l_0^2 + c(\beta)$	l Length of intercept β Angle between test line and surface $c(\beta)$ Function of the planar angle
Surface area for multi-ray designs	$\hat{S} = 2\pi \sum_{j=1}^{2r} l_j^2 \cdot c(\beta)$	l Length of intercept β Angle between test line and surface $c(\beta)$ Function of the planar angle r Number of test lines
Function of the planar angle	$c(\beta) = 1 + \left[\frac{1}{2} \cot \beta \right] \cdot \left[\frac{\pi}{2} - \sin^{-1} \frac{1 - \cot^2 \beta}{1 + \cot^2 \beta} \right]$	β Angle between test line and surface

References

Jensen, E.B., Gundersen, H.J.G. (1987). Stereological estimation of surface area of arbitrary particles. *Acta Stereologica*, 6 (3).