



Neurolucida

Release Notes

Version 2018.1.1 (05/7/2018)

ISSUES RESOLVED

- VSI images: Now loading full dynamic range and correctly reading files that were cropped using Olympus software.
- Zeiss Axioplan 2i or Axioplan 2ie microscope: If the device state toolbar is set to the objective turret, changing the objective to one that requires an oil stop will no longer be a possible cause for crash.
- AxioCam: Updated mcam_zei.h and mcam_zei_ex.h files.
- 3D tree editing: Partial projection now clips trees during editing.