



Version 2025.1.2 (May 2025)

NEW FEATURES AND ENHANCEMENTS

- Improvements made to vessel tracing
- Updates to the Lumina 9 camera settings and monochrome acquires

ISSUES RESOLVED

- Fixed a problem with tracing structures that could cause the software to crash

Version 2025.1.1 (April 2025)

NEW FEATURES AND ENHANCEMENTS

- Support added for the new Lumina 9 camera
- Improvements made to image acquisition with the Hamamatsu camera

ISSUES RESOLVED

- Fixed an issue that could occur when using the Apotome during a slide scan image acquisition
- The storage estimate needed for stitching during a slide scan was fixed

Version 2024.1.3 (November 2024)

LICENSING UPDATE

- Update to offline software licensing. NeuroLucida software requires reactivation on systems that do not have access to the internet. [Instructions for offline software activation](#) are included in the NeuroLucida user guide.

ISSUES RESOLVED

- Fixed a problem with Live Image not working when the program is started
- Fixed an issue with the Vessel Tracing button on the Trace ribbon
- Fixed an issue with the Live Zoom button on the Acquire ribbon