



Neurolucida® 360 software

Release Notes

Version 2020.3.2 (November 2, 2020)

NEW FEATURES AND ENHANCEMENTS

- Added ability to select all dendrites, apical dendrites or axons in 3D environment
- Reporting the thickness of a specific point in the 3D environment
- Added progress bar for Clean-Up
- Added support for additional OME-TIFF format

Neurolucida 360 3D environment

- Fixed preference section for Tree Classification color sets link
- When using subvolume, Puncta detection will be limited to that subvolume
- Fixed the dendrite disappearing after changing the ending type