



NeuroLucida 360

Release Notes

Version 2020.1.1 (12/6/2019)

NEW FEATURES AND ENHANCEMENTS

- Image montage includes new method for 2D image alignment.
- Trace>Select objects: Added 'Select everything of this contour type' option to the right-click menu.
- SPARC: Add help access, improved statuses, icon, and login to Blackfynn access window.
- SPARC Vocabulary Services now launch when a data file or image file is opened, rather than when the software is started.
- SPARC network licensing.
- Large Volume Reconstruction mode enhanced to adjust status warning.

ISSUES RESOLVED

NeuroLucida 360 3D environment

- Fixed inability to use Select All on varicosities with Show Current Section.
- Fixed display of branches in points mode, and varicosities with Show Current Section in some data files.
- Fixed crash when performing right click operations during contour editing.
- Fixed Pan to window center when detecting varicosities on nearest branch.

NeuroLucida 360 overall

- Fixed inability to embed a scalebar.
- Fixed issue when changing contour type from circle contour .
- Improved Batch Tracing user interface.