



NeuroLucida Explorer

Release Notes

Version 2020.2.1 (06/09/2020)

NEW FEATURES AND ENHANCEMENTS

- Added g-ratio for tube wall results in Contour analysis
- Added Proximity analysis of axon to dendrite to the Synapse Marker Detail analysis
- Added report for each individual marker to contour distance
- Added Convex Hull analysis to the Batch analysis

ISSUES RESOLVED

- Fixed issue with Varicosity analysis
- Spine volume now reported to 6 significant digits
- Fixed issue when opening additional data files