



Stereo Investigator

Release Notes

Version 2018.1.1 (05/07/2018)

ISSUES RESOLVED

- VSI images: Now loading full dynamic range and correctly reading files that were cropped using Olympus software.
- Area Fraction Fractionator: The markers' Z placement within a contour now matches the Z value of the section that the contour belongs to.
- SRS Image Stack Series: Preview of SRS grid now available to estimate grid size or display changes.
- Zeiss Axioplan 2i or Axioplan 2ie microscope: If the device state toolbar is set to the objective turret, changing the objective to one that requires an oil stop will no longer be a possible cause for crash.
- AxioCam: Updated mcam_zei.h and mcam_zei_ex.h files.