



# Stereo Investigator®

## Release Notes

### Version 2025.1.2 (May 2025)

#### NEW FEATURES AND ENHANCEMENTS

- Updates to the Lumina 9 camera settings and monochrome image acquisitions

#### ISSUES RESOLVED

- Fixed a problem with drawing contours that could cause the software to crash
- For NL, this will read: Fixed a problem with tracing structures that could cause the software to crash

### Version 2025.1.1 (April 2025)

#### NEW FEATURES AND ENHANCEMENTS

- Support added for the new Lumina 9 camera
- Improvements made to image acquisition with the Hamamatsu camera

#### ISSUES RESOLVED

- Fixed an issue that could occur when using the Apotome during a slide scan image acquisition
- The storage estimate needed for stitching during a slide scan was fixed

### Version 2024.1.3 (November 2024)

#### LICENSING UPDATE

- Update to offline software licensing. Stereo Investigator software requires reactivation on systems that do not have access to the internet. [Instructions for offline software activation](#) are included in the Stereo Investigator user guide.

#### ISSUES RESOLVED

- Fixed a problem with Live Image not working when the program is started
- Fixed an issue with the Live Zoom button on the Acquire ribbon