



Stereo Investigator Cleared Tissue Edition

Release Notes

Version 2020.1.1 (2/20/2020)

NEW FEATURES AND ENHANCEMENTS

- Added controls in the export tracing dialog for exporting images at a specific resolution
- Added an edit option to 'Select everything of this contour type'
- Added a new, fast cursor option in preferences which addresses potential cursor delay issues
- New montage option to 'Order and align using image name'

ISSUES RESOLVED

- Fixed image export issue that exported the background image at a low resolution.
- Fixed crash that occurred when cancelling the tool panel "Close All" button with Image Adjustment window
- Fixed issue with transparency slider in preferences
- Fixed issue when loading RGB nd2 images
- Fixed issue where scaling is over-ridden when images are opened in grid layout
- Fixed issue that 'Show current section' was not respected in the 3D visualization window
- Fixed crash that occurred when changing the tracing options with the circle contour