



Stereo Investigator® Cleared Tissue Edition

Release Notes

VERSION 2023.1.1

NEW FEATURES AND ENHANCEMENTS

- Added ability to create a single section in the serial section manager for support of counting cut slabs of tissue
- Improvements in parameter reporting for Image Volume Fractionator and Image Volume Spaceballs
- Pipelines Ribbon organizes automatic tools in one place
- Batch pipeline for image filtering with minimal effort

ISSUES RESOLVED

- Fixed issue where the information in step 1 of the workflows was not saved to the data file
- Fixed crash during Image Volume SpaceBalls workflow