



Stereo Investigator® Whole Slide Edition

Release Notes

VERSION 2023.1.1

NEW FEATURES AND ENHANCEMENTS

- Updated image processing and image-handling engines that greatly increase the speed of loading, viewing, and saving images
 - Image and data files load noticeably faster than in previous releases.
 - Image adjustments are immediately displayed in both the 2D and 3D windows
- Easier software authorization: the software autofills the online Authorization Request Form with information about your computer system. Just fill in a few fields with information about you and your lab to submit the form.
- Outline Sections function automatically draws contours around tissue sections
- Pipelines Ribbon organizes automatic tools in one place
- Batch pipeline for image filtering multiple files with minimal effort
- Research Resource Identifier (RRID) is now recorded in every xml data file.
- New Freehand Selection tools provide more options for selecting markers or contours.

ISSUES RESOLVED

- Account for missing section in the Cavalieri Estimator is fixed
- Export tracing as an image in monochrome is fixed