



# Vesselucida 360

## Release Notes

### Version 2019.1.2 (7/30/2019)

#### NEW FEATURES AND ENHANCEMENTS

##### Vesselucida 360 overall

- File conversion: Batch converter added (DAT to XML) – to access, type “Convert” in the Quick Launch box.

#### ISSUES RESOLVED

##### Vesselucida 360 3D environment

- Vessels: Automatic tracing creates incorrect loops with the voxel scooping method.
- Contours: Can’t adjust thickness by typing a number; when using the slider to change the thickness then canceling the change, the Cancel operation doesn’t work; contour color may change when selecting a different option under “Show as.”
- Movie: Issues with Playback functionality.

##### Vesselucida 360 overall

- Images: Associated image doesn’t open with data file after adjusting the histogram but not saving the changes.
- Image adjustment: Can’t scroll down to view all channels.
- Contours: Movements can be inconsistent.
- AutoMove: Size of AutoMove box changes when Auto-contour is selected.