



Vesselucida 360

Release Notes

Version 2019.1.3 (10/18/2019)

NEW FEATURES AND ENHANCEMENTS

- SPARC Vocabulary Services: Can use symbols in Subject ID field.
- Trace>Select objects: Added 'Select everything of this contour type' option to the right-click menu.

ISSUES RESOLVED

Vesselucida 360 3D environment

- Image: Image no longer visible when using Image Slice mode.

Vesselucida 360 overall

- Dockable windows: 'Close all' button doesn't close Image Adjustment window and may cause crash.
- SPARC Vocabulary Services: Search Contours bar doesn't let the user search for the first contour; list from previous session is displayed and FMA Id is not written to file.
- Trace>Contour selection: "ı »ç" characters are appended to the beginning of the first contour name upon uploading a .csv file; contour names hidden from Contour Settings window still displayed in Contour Selection window.
- SPARC Blackfynn Portal: Issues with configuring Blackfynn agent; after selecting 'Remember me' option when signing in for the first time, cannot sign in again.