



Vesselucida 360

Release Notes

Version 2020.1.1 (02/06/2020)

NEW FEATURES AND ENHANCEMENTS

- Added a new method for 2D image alignment in Image montage.
- Added 'Select everything of this contour type' option to the right-click menu in Trace>Select objects.
- Enhanced Large Volume Reconstruction mode to adjust status warning.

SPARC USERS

- Added help access and improved the Blackfynn access window statuses, icon, and login process.
- can you say something like, "Improved the software start speed by launching SPARC Vocabulary Services when a data file or image file is opened, rather when the software is started.
- Made SPARC Network Licensing available.

ISSUES RESOLVED

Vesselucida 360 3D environment

- Fixed display of vessels in points mode with Show Current Section in some data files.
- Fixed crash that occurred when using right-click operations during contour editing.

Vesselucida 360 overall

- Fixed inability to embed a scalebar.
- Fixed issue when changing contour type from circle contour.