



WormLab[®]

Release Notes

VERSION 2022.1.1 JANUARY 5, 2022

HOW TO DOWNLOAD

If you are an existing WormLab customer, you can download the latest WormLab software from the MBF Download center (<http://support.mbfbioscience.com>).

Microsoft Windows 32 bit, 64 bit and Mac OSX versions are available for WormLab 2022.

On Windows, if you have a 64 bit machine, we highly recommend using the 64 bit version of WormLab.

NEW FEATURES AND ENHANCEMENTS

Video Recording

- New, more efficient video encoder improves encoding speed and decreases file size
- Increased available buffer to support fast frame-rate recording using high resolution cameras
- Added support for latest Basler cameras. Pylon version 6.2.0 or later is required.

General

- Added support for macOS Big Sur
- Added option to specify a fixed com port for systems where auto scanning ports may fail
- Measuring tool for scaling may be drawn at any angle

ISSUES RESOLVED

- Fixed an issue in the curvature map analysis where incorrect curvature values were displayed in the data table in degree mode
- Fixed an issue in the computation of the standard deviation in Group Analysis